



Deta  
Elis  
Holding  
Health  
Innovation  
Human

A wellness device for  
enhancing our own protective  
forces and correcting  
the natural electromagnetic  
radiation in natural bio system

DeVita AP+  
Model Mini

Methodological recommendations



## Introduction

A unique device which has no analogues in the world. Each person needs to maintain a high rate of health in the modern world. The device is used to relieve the stressful burden caused by the presence of parasitic forms in the human body: viruses, bacteria, fungi, protozoa, helminths. And also the presence of endo- and exo- toxins that block the work of the main excretory systems, which is one of the main causes of reducing the state of health.

## Method of work

Its action is based on the results of large-scale research by the scientist and engineer R. Rife and his supporters. Scientists have proved that every living form has its own electromagnetic radiation. Waves of a certain frequency, reproduced by the device, resonate with the own radiation of the microorganisms, enhancing their natural vibrations. Intensive fluctuations become disastrous for parasites. The complex AP+ Mini device has separate vibration frequencies for certain types of viruses, bacteria, helminths, protozoa and fungi, and, in addition, allows the deep purification of the body from toxins.

## Features of the device

The device includes 27 (7 of them new) complexes consisting of manual and automatic programs. Each user can create their own unique complexes (the device can hold up to 50 in total) and has the functional to download new complexes and programs that will be created by the company's specialists.

## Time of use

The DeVita AP+ Mini device should be used in the afternoon or at night, placing it under the pillow with the backside facing yourself.

To maintain a high wellness factor, to improve the state of health and for best adaptation, we recommend using the DeVita Energy 11 Mini device in the morning and the DeVita Ritm+ Mini device during the day.

These Methodological Recommendations contain a list of all programs included in the following complexes, application diagrams of the DeVita AP+ Mini device for complexes No. 4 - 20, and a description of automatic programs. Detailed application schemes of the devices are given in the "Complex Approach" catalog. A description of all manual programs and their application is given in the brochure "List of manual programs of the DeVita AP + Mini".

## Complexes

- |                               |                      |                      |
|-------------------------------|----------------------|----------------------|
| 1. General list of programs 1 | 10. Sport performers | 19. Healthy teeth    |
| 2. General list of programs 2 | 11. Clean skin       | 20. Stop Coronavirus |
| 3. General list of programs 3 | 12. Stop parasites   | 21. First aid        |
| 4. Heart and circulation      | 13. Stop allergy     | 22. Gastritis +      |
| 5. Digestion without problems | 14. Men's health     | 23. Healthy skin     |
| 6. Stop stress                | 15. Women's health   | 24. Healthy back     |
| 7. Stop insomnia              | 16. Slim figure      | 25. Stop infection   |
| 8. No cold                    | 17. Healthy Cell     | 26. Health 60+       |
| 9. Movement without pain      | 18. Lyme Disease     | 27. Healthy child    |

### 1. Complex "General list of programs 1"

#### Manual programs

1. Anti-Cold
2. Anti-Flu
3. Antivirus
4. Healthy gastrointestinal
5. Intimate health
6. Stop candida
7. Against flu - General
8. No herpes
9. Without toxins
10. No herpes simplex
11. No herpes zoster
12. No mould & fungi
13. No chlamydia
14. No parasites, Digestive tract
15. No parasites, General
16. No protozoans
17. No staphylococcus
18. No streptococcus
19. Weight Control
20. Against Lyme Disease 1
21. Against Lyme Disease 2
22. Antivirus Gastrointestinal
23. Healthy Cell 1
24. Healthy Cell 2
25. Stop Amoeba
26. Stop Penicillium
27. No Trichomonas
28. Drainage - General

29. Healthy lymph

#### Automatic programs

30. Mode 1-automatic
31. Mode 2-automatic
32. Mode 3-automatic
33. Mode 4-automatic
34. Cleansing-automatic
35. No parasites-automatic
36. No flu-automatic
37. Stop herpes S-automatic
38. Stop herpes Z-automatic
39. Without cold-automatic
40. Support immune system-Automatic
41. Weight Reduction-Automatic
42. Reduce Points-automatic
43. Cell Protection 1 - Automatic
44. Cell Protection 2 - Automatic
45. Cell Protection 3 - Automatic
46. Cell Protection 4 - Automatic
47. Lyme 1 -Automatic
48. Lyme 2 - Automatic
49. Lyme 3 - Automatin
50. Lyme 4 - Automatic

## 2. Complex “General list of programs 2”

### Manual programs

1. Without Salmonella
2. Without Yersinia
3. Without Proteus
4. Without Klebsiella
5. Without Campylobacter
6. Without Clostridia
7. Without Shigella
8. Without Rotavirus
9. Without Trematodes – General 1
10. Without trematodes – General 2
11. Without Enterococcus
12. Without nematodes and ascarididae  
- New
13. Stop adenovirus - General
14. Stop Coronavirus
15. Stop Rickettsia
16. Without Papilloma - General
17. Without Papilloma
18. Without Leptospira
19. Without Cryptococcus
20. Without Toxoplasma
21. Without Ureaplasma
22. Without Echinococcus
23. Without Pseudomonas aeruginosa
24. Help against poisoning
25. No bacteria
26. No fungi
27. No helicobacter
28. No hepatitis A
29. No schistosomes (blood-flukes) & trematods
30. No nematodes and ascarides

## 3. Complex “General list of programs 3”

### Manual programs

1. Healthy Gums
2. Mucous membrane Regeneration
3. Healthy mouth
4. Without Fistula
5. Stop Injury (teeth)
6. CNS Without Problems
7. Healthy Nerves
8. Clean Sinuses - General
9. Clean Sinuses
10. Ears Without Inflammation
11. Throat Without Pain
12. Against Cold - General
13. Anti-inflammation - General
14. Stop Migraine
15. Antiseptic
16. Healthy Bile
17. Drainage Antitox
18. Mexican Detoxication
19. No allergies

### Automatic programs

20. Restoration 1 - Automatic
21. Restoration 2 - Automatic
22. Healthy Gums - Automatic
23. Nerve Regeneration - Automatic

**Attention!** If during the use of the programs of the first month nausea, headache, perspiration, or other symptoms of intoxication are observed, then additional application of this or other mode has to be carried out within the next 3 months. Each of them is involved in the course taking into account the life cycle of the infectious organism.

So:

- Modes 1-3: apply 14 days;
- Mode 4: 30 days.

The device is used daily at night

## 4. Complex “Heart and circulation”

Applied for problems with the cardiovascular and lymphatic systems.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1- automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month:

- 1<sup>st</sup> -2<sup>nd</sup> week — program “No schistosomes (blood-flukes) & trematods” daily use in the evening after the programs of DeVita Ritm Mini.
- 3<sup>rd</sup> -4<sup>th</sup> week — program “No chlamydia” daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — “Cleansing” daily at night.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: “No parasites-automatic” daily at night time for 1 week.

## 5. Complex “Digestion without problems”

For chronic problems with the organs of the gastrointestinal tract and as an additional tool in case of poisoning.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1- automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month:

- 1<sup>st</sup> -2<sup>nd</sup> week — program “No schistosomes (blood-flukes) & trematods” daily use in the evening after the programs of DeVita Ritm Mini.
- 3<sup>rd</sup> -4<sup>th</sup> week — program “Drainage antitox” daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — “Cleansing” daily at night.

#### If necessary:

“No helicobacter” — daily in the evening for a month.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: “No parasites-automatic” daily at night time for 1 week.

#### In case of poisoning:

- “Mode 3 - automatic”
  - program “Help against poisoning”
- Daily for 10 days, alternating with each other.

## 6. Complex “Stop stress”

In case of stress, fatigue, reduced performance, depressive states.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1 - automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month: “Cleansing” daily at night.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: “No parasites-automatic” daily at night time for 1 week.

## 7. Complex “Stop insomnia”

In case of a state of nervous excitement, trouble falling asleep, sleep disorders.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1 - automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month: "Cleansing" daily at night.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: "No parasites-automatic" daily at night time for 1 week.

## 8. Complex "Active immune system"

For the prevention of colds, as well as relief from bacterial and viral infections.

In the season of colds, in case of hypothermia:

- "Support immunity –automatic" 1-2 times a week

At the first symptoms of the common cold diseases (other than flu) during the first two days:

- "Without cold – automatic"

Then within 5 to 7 days alternating with each other:

- "Mode 2 - automatic"

- program "Anti-Cold"

At the first signs(symptoms) of the flu, for the first 2 days alternating with each other:

- program "Antivirus"

- "No flu –automatic"

Further, daily for 5 days: "No flu – automatic"

Upon completion, use "Mode 2 - automatic" daily for 5-7 days, also if necessary the program "No bacteria".

## 9. Complex "Movement without pain"

In case of problems with the musculoskeletal system: pain, inflammation of the joints, muscles, spine or injuries.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — "Mode 1 - automatic"

- 2<sup>nd</sup> week — "Mode 2 - automatic"

- 3<sup>rd</sup> week — "Mode 3 - automatic"

- 4<sup>th</sup> week — "Mode 4 - automatic"

2<sup>nd</sup> month:

- "Cleansing" daily at night.

- program "No bacteria" and "Drainage antitox". Daily in the evening after the DeVita Ritm Mini programs for 14 days.



3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: "No parasites-automatic" daily at night time for 1 week.

## 10. Complex "Sport performers"

With intense physical activity, injuries, as well as to reduce the toxic load and a fast recovery after workouts.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — "Mode 1 - automatic"
- 2<sup>nd</sup> week — "Mode 2 - automatic"
- 3<sup>rd</sup> week — "Mode 3 - automatic"
- 4<sup>th</sup> week — "Mode 4 - automatic"

2<sup>nd</sup> month: "Cleansing" daily at night.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: "No parasites-automatic" daily at night time for 1 week

## 11. Complex "Clean skin"

With inflammatory processes on the skin, fungal burdening of hair and nails, herpetic weights (shingles).

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — "Mode 1 - automatic"
- 2<sup>nd</sup> week — "Mode 2 - automatic"
- 3<sup>rd</sup> week — "Mode 3 - automatic"
- 4<sup>th</sup> week — "Mode 4 - automatic"

2<sup>nd</sup> month: "Cleansing" daily at night.

### If necessary:

- in case of fungal burdening program "No fungi" daily in the evening after the programs of DeVita Ritm Mini

3<sup>rd</sup> month:

- necessary use of the programs of the first month
- in case of herpetic burdening: "Stop herpes S-automatic" or "Stop herpes Z-automatic" at night time.

4<sup>th</sup> month: "No parasites-automatic" daily at night time for 1 week.

## 12. Complex “Stop parasites”

To reduce health risks of: pet owners, gardeners, travelers, families with children under 12 years. In case of appetite changes, unexplained pain, inflammation, appearance of allergies.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1 - automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month: “Cleansing” daily at night.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: “No parasites-automatic” and “Without helminths – automatic”, alternating with each other every day at night for 1 month.

Additionally:

In case of cough, program “No bacteria” — daily in the evening after the programs of DeVita Ritm Mini, for 14 days.

In case of herpetic weights:

- program “No herpes simplex” at the first sign of a rash every 3 hours, then every day 1-2 times a day “Stop herpes S-automatic” at night for 1 month.
- program “No herpes zoster” at the first sign of a rash every 3 hours, then every day 1-2 times a day “Stop herpes Z-automatic” at night for 1 month.

When visiting countries in Asia and Africa, to prevent viral load of the liver, use the program – “No hepatitis A” – daily 1-2 times a day for a month.

## 13. Complex “Stop allergy”

With various allergic reactions to: food, dust, pets, pollen of plants.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1 - automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

Additionally: program “No allergies” in the evening after the programs of DeVita Ritm Mini.

2<sup>nd</sup> month:

- 1<sup>st</sup> -2<sup>nd</sup> week — program “No schistosomes (blood-flukes) & trematods” daily use in the evening after the programs of DeVita Ritm Mini.
- 3<sup>rd</sup> -4<sup>th</sup> week — program “No allergies” daily use in the evening after the programs of DeVita Ritm Mini.
- 1-4 week — “Cleansing” daily at night.

3<sup>rd</sup> month:

- necessary use of the programs of the first month
- in case of herpetic burdening: “No allergies” daily at night.

4<sup>th</sup> month:

- “No parasites-automatic” daily at night time for 1 week.
- program “No allergies” in the evening after the programs of DeVita Ritm Mini.

**Attention!** In case of pollinosis, the procedure should be started 2 months before the flowering of plants.

A stable positive result for solving problems with allergies can be obtained with regular use of the programs of the complex for 6 months.

## 14. Complex “Men’s health”

In case of psychological and functional problems with the organs of the male urogenital system.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — “Mode 1 - automatic”
- 2<sup>nd</sup> week — “Mode 2 - automatic”
- 3<sup>rd</sup> week — “Mode 3 - automatic”
- 4<sup>th</sup> week — “Mode 4 - automatic”

2<sup>nd</sup> month:

- “Cleansing” daily at night.
- program “No chlamydia” in the evening after the programs of DeVita Ritm Mini.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: “No parasites-automatic” daily at night time for 1 week

## 15. Complex “Women’s health”

In case of functional disorders of the female urogenital system, for harmonizing the work of the organism in premenstrual and menopausal period.

1<sup>st</sup> month of daily use, during night time, according to the following scheme:

- 1<sup>st</sup> week — "Mode 1 - automatic"
- 2<sup>nd</sup> week — "Mode 2 - automatic"
- 3<sup>rd</sup> week — "Mode 3 - automatic"
- 4<sup>th</sup> week — "Mode 4 - automatic"

2<sup>nd</sup> month:

- "Cleansing" daily at night.
- program "No chlamydia" in the evening after the programs of DeVita Ritm Mini.

3<sup>rd</sup> month: necessary use of the programs of the first month

4<sup>th</sup> month: "No parasites-automatic" daily at night time for 1 week

## 16. Complex "Body figure"

For soft and comfortable weight loss, regulation of appetite.

**Attention!** To obtain the most effective result before using this complex, we recommend first to reduce the possible parasitic and toxic load, using the basic programs ("Mode 1-4 automatic", "Cleansing", "No parasites-automatic") of the complex "Stop parasites" within 3-4 months.

4<sup>th</sup>-6<sup>th</sup> month: programs "Weight Reduction-automatic" and "Reduce Points – automatic" every day, alternating with each other at night

## 17. Complex "Healthy Cell"

It is used to prevent cell degeneration and the formation of pathological neoplasms.

**Attention!** This complex is recommended to be used after a 2-month basic cleansing according to the regimen of the "Digestion without problems" complex.

3<sup>rd</sup> month daily at night as follows:

- 1<sup>st</sup>-2<sup>nd</sup> week — "Cell Protection" and "Cell Protection 2" alternate every other day. The program start time is not later than 17:00.
- 3<sup>rd</sup> -4<sup>th</sup> week — "Cell Protection 3" and "Cell Protection 4" alternate every other day. The program start time is not later than 17:00.

4<sup>th</sup> month:

- 1<sup>st</sup> day — program "Cleansing-automatic" at night
- 2<sup>nd</sup> day — program "Without papilloma - General" at night or at the evening.

**Programs of 3<sup>rd</sup> and 4<sup>th</sup> months are to repeat for another 2 cycles!**

## 18. Complex “Lyme Disease”

It is used to maintain a good health coefficient after suffering tick-borne borreliosis (Lyme disease).

**Attention!** This complex is recommended to be used after a 2-month basic cleansing according to the regimen of the “Digestion without problems” complex.

**3<sup>rd</sup> month** daily at night as follows:

- 1<sup>st</sup>-2<sup>nd</sup> week — “Lyme 1 - Automatic” and “Lyme 2 - Automatic” alternate every other day. The program start time is not later than 17:00.
- 3<sup>rd</sup>-4<sup>th</sup> week — “Lyme 3 - Automatic” and “Lyme 4 - Automatic” alternate every other day. The program start time is not later than 17:00.

**4<sup>th</sup> month:**

- 1<sup>st</sup> day — program “Cleansing-automati” at night.
- 2<sup>nd</sup> day — program “Nerve Regeneration - Automatic” at night or at the evening

**Programs of 3<sup>rd</sup> and 4<sup>th</sup> months are to repeat for another 2 cycles!**

## 19. Complex “Healthy teeth”

It is used after dental interventions and in chronic problems of the oral cavity.

**Scheme 1.** Use after dental interventions (removal, implantation, treatment, etc.). Duration of use is 2 weeks:

- 1<sup>st</sup> day — program “Restoration 1 Automatic” in the evening after the DeVita Ritm Mini programs. “Mode 2-automatic” at night.
- 2<sup>nd</sup> day — program “Restoration 2 Automatic” in the evening or at night.

In the presence of abscesses, fistulas add the program daily “Without Fistula” daily.

**Scheme 2.** Use for chronic oral problems:

**One month** every day at night according to the regimen of the “Digestion without problems” complex.

**2<sup>nd</sup> month:**

- 1<sup>st</sup> day — program “Healthy Gums Automatic” at night.
- 2<sup>nd</sup> day — program “Cleansing-automati” at night.

**3<sup>rd</sup> month** if needed: 1<sup>st</sup> month programs.

**4<sup>th</sup> month:**

- “No parasites-automatic” daily at night time for 1 week.

## 20. Complex “Stop Coronavirus”

To prevent coronavirus infection, as well as alleviate the condition with bacterial and viral loads.

In order to reduce the risk of infection 2-3 times a week:

- “Without cold-automatic” at night.
- “Stop coronavirus new” after DeVita Ritm programs.

At the first signs of the disease, alternating daily with each other every other day:

During the 1<sup>st</sup> week and at night:

- “Without cold-automatic”
- “Healthy breathing automatic”

Additionally:

- “Restoration 2 - Automatic” after the DeVita Ritm program “Healthy day automatic”.

During the 2<sup>nd</sup> week at night:

- “Base 2 automatic”
- “Healthy breathing automatic”

Additionally after the DeVita Ritm program “Health”:

- “Without cough”
- “Without pneumonia”

To maintain health after infection, daily alternating with each other every other day for 2 weeks at night:

- “Restoration 2 - Automatic”
- “Base 4 automatic”

## 21. Complex “First aid”

The effect of the programs of the complex is aimed to quickly eliminate unpleasant sensations: nausea, pain, inflammation.

### Manual programs

1. Low back without inflammation
2. Stop pain (lower back)
3. Stop pain (injury)
4. Stop pelvic inflammation
5. Stop spasm (muscles)
6. Without motion sickness
7. Stop spasm (Intestines)
8. Restoration of blood circulation
9. Healing
10. Male power
11. Without allergies new
12. Stop conjunctivitis

13. Strong nerves
14. Good sleep
15. Without hemorrhoids
16. Head circulation
17. For women
18. Stop fatigue

### Automatic programs

19. Without motion sickness -automatic
20. Sos! Allergy - Automatic
21. Restoration 3 -automatic
22. Deep sleep-automatic
23. Stop fatigue-automatic
24. Anti-inflam. (female) -automatic
25. Anti-inflam. (male) -automatic
26. Low back without pain -automatic

## 22. Complex “Gastritis +”

The effect of the programs of the complex is aimed at improving the work of the GI organs: It reduces inflammation, restores mucous membranes, eliminates discomfort and ensures normal functioning of the stomach, all parts of the intestines, pancreas and liver.

### Manual programs

1. Without fluke
2. Calm intestines
3. Healthy Stomach
4. Duodenitis, stop!
5. Healthy pancreas
6. Blood sugar under control
7. Stop diabetes
8. Small intestines
9. Intestines without disorders
10. Without constipation
11. Without motion sickness
12. Stop spasm (Intestines)
13. Healthy spleen

### Automatic programs

14. Without motion sickness -automatic
15. Active Gastro -automatic
16. Bowel without irritation
17. Without constipation -automatic

## 23. Complex “Healthy skin”

The effect of the programs of the complex is aimed at fighting pathogenic bacteria and fungi of the skin, reducing inflammation, increasing skin immunity, improving the external condition of the skin, hair, and nails.

### Manual programs

1. Stop acne bacteria
2. Without microsporium
3. No onychomycosis
4. Zygomycetes, stop!
5. Nails without fungus
6. Without trichophyton

7. Mycosis fungoides
8. Without epidermophyton
9. Without actinomycetes
10. Fungus-free skin
11. Solar keratosis
12. Stop warts

13. Stop pimples
14. Fingers without inflammation
15. No scleroderma
16. Erysipelas, stop
17. Relieve psoriasis
18. Leprosy
19. Skin without acne
20. Without toxins

### Automatic programs

21. Skin without pimples -automatic
22. Healthy nails -automatic
23. Stop skin fungus -automatic

## 24. Complex “Healthy back”

The effect of the complex programs is aimed at improving the condition of the spine and back muscles, reducing pain sensations, reducing inflammation and swelling. It facilitates the ease of movement.

### Manual programs

1. Low back without inflammation
2. Straight spine
3. Stop pain (low back)
4. Stop spasm (muscles)
5. Muscles without tension

### Automatic programs

6. Low back without pain -automatic
7. Healthy back -automatic

## 25. Complex “Stop infection”

The effect of the programs of the complex is aimed at destroying pathogenic viruses and bacteria that are the pathogenic agents of the most common infections, as well as the removal of toxins from the body.

### Manual programs

1. Stop chickenpox
2. Stop rubella
3. Without mumps
4. Stop scarlet fever
5. Stop measles
6. No diphtheria
7. Stop pertussis
8. Stop gardnerella
9. Stop gonococcus
10. Stop parvovirus
11. Without tuberculosis
12. Without brucella
13. Stop cholera
14. Without toxins



## 26. Complex “Health 60+”

The effect of the programs of the complex is aimed at maintaining the health of elderly people, taking into account the main problems that are typical for this category of the population. It will help in achieving an active longevity.

### Manual programs

1. Stop dementia
2. Stop Parkinson
3. Head circulation
4. Strong nerves
5. Good sleep
6. Solar keratosis
7. Stop fatigue
8. Without constipation
9. Anti-Cold
10. Antivirus
11. Healthy gastrointestinal
12. Low back without inflammation
13. Stop pain (lower back)
14. Stop pain (injury)
15. Stop pelvic inflammation
16. Without allergies new
17. Stop conjunctivitis
18. Strong nerves
19. Good sleep
20. Without hemorrhoids
21. Calm intestines
22. Healthy Stomach
23. Duodenitis, stop!
24. Healthy pancreas
25. Blood sugar under control
26. Anti-inflammation-General
27. Stop Migraine
28. Antiseptic
29. Healthy Bile
30. Stop spasm (muscles)
31. Restoration of blood circulation
32. Healing
33. Healthy Gums
34. No onychomycosis
35. Without trichophyton
36. Nails without fungus
37. Without microsporum
38. Bacterial pneumonia

39. Stop coronavirus NEW
40. Anti-virus Gastrointestinal
41. Without toxins
42. Healthy lymph

### Automatic programs

43. Stop brain aging -automatic
44. Deep sleep-automatic
45. Healthy nails -automatic
46. Low back without pain -automatic
47. Active Gastro -automatic
48. Without cold-Automatic
49. Healthy breathing-Automatic
50. Restoration 3 -automatic

## 27. Complex “Healthy child”

The effect of the programs of the complex is aimed at maintaining the health of children from 0 to 12 years of age, taking into account the age particularities. It will help to strengthen immunity, better adapt to new conditions (kindergarten, school).

### Manual programs

1. Stop chickenpox
2. Stop rubella
3. Without mumps
4. Stop scarlet fever
5. Stop measles
6. No diphtheria
7. Stop pertussis
8. No hepatitis A
9. Anti-Cold
10. Antivirus
11. Help against poisoning
12. Healthy gastrointestinal
13. Bacterial pneumonia
14. No streptococcus
15. No staphylococcus
16. No parasites, General
17. Without nematodes and ascarididae-New
18. No protists
19. Without toxins
20. Anti-inflammation-General
21. Antiseptic
22. Without allergies new
23. Without motion sickness
24. Stop spasm (Intestines)
25. Healthy Stomach
26. Duodenitis, stop!
27. Without constipation
28. Healthy Bile
29. Intestines without disorders
30. Straight spine
31. Stop spasm (muscles)
32. Muscles without tension

33. Against Cold-General
34. Ears Without Inflammation
35. Throat Without Pain
36. Against flu-General
37. Anti-Flu
38. Stop pain (injury)
39. Without cough
40. Clean Sinuses-General

### Automatic programs

41. Without cold-Automatic
42. Mode 2-automatic
43. Without motion sickness -automatic
44. Healthy back -automatic
45. Restoration 2-Automatic
46. Without constipation -automatic

With regular, correct use of DeVita Mini devices any user will receive a unidirectional positive result of effective recovery, youth and longevity.

## Description of automatic programs

Name	Effect of the program	Programs	Working time
Mode 1 Automatic	It inhibits the activity of viruses, inhibits the vital activity of protozoan microorganisms. Optimal disease prevention when changing climatic zones, visiting Asian countries.	Antivirus, No protozoans, Without toxins	6:47
Mode 2 Automatic	Relieves toxic load caused by streptococcal and staphylococci. It is effective both for colds and for the prevention of coccal infections.	No staphylococcus, No streptococcus, Without toxins	7:53
Mode 3 automatic	It contributes to the removal of bacterial aggravation of the digestive tract and the genitourinary system.	Healthy gastrointestinal, Intimate health, Without toxins	8:26
Mode 4 Automatic	It removes mycotoxins formed by mold and yeast. Reduces the risk of accumulation of carcinogenic factors	Stop candida, No mould & fungi, Without toxins	7:47
Cleansing-automatic	It promotes the elimination of toxins, the normalization of lymph outflow, and the strengthening of lymphatic vessels. Helps the removal of excess fluid from the intercellular space, has an antiseptic effect.	Healthy lymph, Drainage - General	7:10
No parasites Automatic	Reduces the stress parasitic burden of the liver, pancreas, intestines, as well as the toxic load caused by the vital activity of these forms.	No parasites, Digestive tract, No parasites, General, Without toxins	7:38
No flu - automatic	It suppresses the activity of various strains of the influenza virus, helps relieve symptoms of intoxication and alleviates the human condition during the illness.	Anti-Flu, Against flu - General, Without toxins	6:38

Stop herpes S-automatic	It prevents the development of herpes simplex and its toxic effect on the body: the formation of vesicles on the mucous membranes, skin, reduces itching.	No herpes, No herpes simplex, Without toxins	7:35
No herpes Z-automatic	It reduces the activity of zoster herpes viruses, reduces the toxic signs of its manifestation on the surface of the skin and pain.	No herpes zoster, No herpes, Without toxins	5:47
Without cold automatic	Prevents the development of respiratory viruses, reduces symptoms of intoxication, leads to a decrease in inflammatory reactions and the likelihood of complications	Antivirus, Anti-Cold, Without toxins	8:09
Support immune system-Automatic	It blocks the activity of various viruses, their distribution in the human body, promotes the activation of immunity. It is effective as prevention of colds infections during hypothermia and during the period of mass diseases.	Antivirus, Without toxins	5:03
Weight Reduction-Automatic	Reduces parasitic load and destroys viruses that are the cause of excess weight. It normalizes the functioning of the autonomic nervous system, maintains the body tone and improves the psychoemotional state.	Weight Control, Without toxins	4:28
Reduce Points-automatic	It cleanses the body of toxins, normalizes lymphatic drainage, supports the immune system and increases vital energy.	Healthy lymph, Without toxins	4:31
Cell Protection 1 -Automatic	The program helps to relieve viral, bacterial and toxic loads, increasing the body's defenses.	Antivirus + No streptococcus + No herpes simplex + Without toxins	12:09
Cell Protection 2 – Automatic	Prevents the effects of carcinogenic factors on the human body, thereby reducing the risk of cell degeneration.	Healthy Cell 1 + Healthy Cell 2 +Antivirus Gastrointestinal + No Trichomonas	13:54

Cell Protection 3 – Automatic	The program helps to resist viral, bacterial and protozoa infections, reducing the degree of intoxication.	No herpes zoster + No staphylococcus + No chlamydia + Intimate health + Stop Amoeba	12:42
Cell Protection 4 – Automatic	The program reduces the fungal toxic load, increasing your own protection and reduces the likelihood of tumor processes.	Healthy Cell 1 + Healthy Cell 2 + Stop candida + Stop Penicillium	12:48
Lyme 1 - Automatic	The program helps to remove viral, bacterial and toxic loads, increasing the body's defenses.	Antivirus + No streptococcus + No herpes simplex + Without toxins	12:09
Lyme 2 - Automatic	The program works to destroy pathogens that have been ingested by a tick bite and to reduce the total toxic load.	Against Lyme Disease 1 + Against Lyme Disease 2 + Antivirus Gastrointestinal + No Trichomonas	11:39
Lyme 3 - Automatic	The program helps to resist viral, bacterial and protozoa infections, reducing the degree of intoxication.	No herpes zoster + No staphylococcus + No chlamydia + Intimate health + Stop Amoeba	12:42
Lyme 4 - Automatic	The program reduces the fungal toxic load, increasing the body's defenses while fighting borrelia.	Against Lyme Disease 1 + Against Lyme Disease 2 + Stop candida + Stop Penicillium + Without toxins	11:51

Restoration 1 - Automatic	The program is aimed to relieve pain after dental interventions, repair damaged mucosa, reduce swelling, and accelerate tissue regeneration.	Stop Injury (teeth) + Healthy mouth	6:33
Restoration 2 - Automatic	The program reduces the influence of pathogenic factors, reducing the degree of the inflammatory process. Promotes speedy recovery.	Antiseptic+Anti-inflammation - General	6:09
Healthy Gums - Automatic	The program helps to reduce the inflammatory process of the oral cavity, improves trophic tissue, strengthens local immunity.	Healthy Gums+Mucous membrane Regeneration	6:48
Nerve Regeneration - Automatic	The program is aimed at improving the functioning of the nervous system in case of damage by pathogens that have entered the body due to a tick bite. It reduces swelling, helps restore the integrity of the nerve fiber and nerve connections of various parts of the nervous system.	Healthy Nerves+CNS Without Problems	4:51
Healthy Breathing	The program aims to suppress the activity of viruses and bacteria that cause respiratory tract diseases, helps to reduce the likelihood of exacerbations, reduce intoxication, increase the body's defenses.	Stop coronavirus NEW + Bacterial pneumonia + Without toxins	10:37
Without motion sickness -automatic	The program helps relax the intestinal wall, reduces pain, nausea, and promotes better health. <b>Application:</b> use preferably 30 min before travel.	Without motion sickness Stop spasm (Intestines)	4:35
Sos! Allergy - Automatic	The program reduces allergic manifestations, promotes better breathing by relieving spasm, reducing inflammation, relieving swelling and itching.	Without allergies new Stop conjunctivitis	6:14

Sos! Allergy - Automatic	<p><b>Application:</b> daily for the entire period of exposure to allergens or as needed, in conjunction with the corresponding programs on DeVita Ritm (see "Stop Allergy" complex).</p>	Without allergies new Stop conjunctivitis	6:14
Restoration 3 -auto- matic	<p>The program improves blood supply to organs and tissues, activates the regeneration process especially after infectious diseases or injuries.</p> <p><b>Application:</b> daily for a month, if necessary use the program until full recovery, in conjunction with DeVita Ritm programs (see Heart and blood vessels complex).</p>	Restoration of blood circulation Healing	2:53
Deep sleep- automatic	<p>The program normalizes the psychological state, reduces the excitability of the nervous system, improves mood and promotes restful sleep.</p> <p><b>Application:</b> daily for a month, in conjunction with DeVita Ritm programs (see "Stop Stress" complex) or "Stop insomnia"</p>	Strong nerves Good sleep	6:20
Stop brain aging - automatic	<p>The program helps maintain brain function, slows down the process of neuron destruction, helps improve blood supply to the tissues of the head organs, thereby increasing the flow of oxygen and nutrition to them.</p> <p><b>Application:</b> daily for a month, in conjunction with the corresponding programs on program on DeVita Ritm (see Heart and blood vessels complex). Repeat the cycle regularly every quarter.</p>	Stop dementia Stop Parkinson Blood circulation in the head	9:36
Stop fatigue- automatic	<p>The program reduces the agitation of the nervous system, helps to normalize the psychological state, improves blood supply to the tissues of the organs of the head, helps to increase performance.</p>	Strong nerves Blood circulation in the head Stop fatigue	9:09

Stop fatigue-automatic	<b>Application:</b> daily, in conjunction with DeVita Ritm programs (see "Stop stress" complex).	Strong nerves Blood circulation in the head Stop fatigue	9:09
Anti-inflammation (female) -automatic	The program reduces the inflammatory process, relieves pain, helps to restore hormonal balance, normal functioning of the reproductive system. <b>Application:</b> daily for a month, in conjunction with the DeVita AP and DeVita Ritm programs for relevant infectious agents (see "Women's Health" complex).	Stop pelvic inflammation For women	6:35
Anti-inflammation (male) -automatic	The program helps to reduce the inflammatory process, restore the normal functioning of the reproductive system, reduces the chance of complications (impotence, infertility). <b>Application:</b> daily for a month, in conjunction with the DeVita AP and DeVita Ritm programs for relevant infectious agents (see "Women's Health" and "Men's Health" complexes).	Stop pelvic inflammation Male power	5:56
Active Gastro -automatic	The program helps to relieve pain, reduce inflammation, restore mucosa, normalize the stomach, gallbladder, and pancreas. <b>Application:</b> daily for a month, in conjunction with the DeVita Ritm program "Stomach+ automatic" and the programs of the complex "Stop stress" or "Digestion without problems". In order to prevent exacerbations, the cycle is carried out 2-3 times a year.	Healthy Stomach Duodenitis, stop! Healthy pancreas	9:21
Bowel without irritation	The program reduces pain, meteorism, bloating, relieves inflammation of the mucosa, normalizes the motility of the GI tract.	Calm intestines Intestines without disorders	7:56



Bowel without irritation	<p><b>Application:</b> daily for a month, in conjunction with the corresponding programs on DeVita Ritm (see "Digestion without problems" complex).</p>	Calm intestines Intestines without disorders	7:56
Without constipation -automatic	<p>The program helps to relieve spasms, reduce pain, normalizes the motility of the large intestine.</p> <p><b>Application:</b> daily for a for one month, in conjunction with the corresponding programs on DeVita AP and DeVita Ritm (see "Digestion without problems", "Stop stress" complexes).</p>	Without constipation Stop spasm (Intestines)	3:44
Skin without pimples -automatic	<p>The program fights pathogenic microflora, reduces inflammation, removes toxins from the body, regulates the barrier function of the skin, improves its condition.</p> <p><b>Application:</b> daily for a month, in conjunction with the "Mode 2-automatic" program on DeVita AP and the corresponding programs on DeVita Ritm (see Clear Skin complex)..</p> <p>In case of skin predisposition to acne formation, in order to prevent inflammation, it is desirable to use the program every quarter according to the above scheme.</p>	Stop acne bacteria Without toxins Skin without acne Stop pimples	5:01
Healthy nails - automatic	<p>The program causes the death of pathogenic fungi, reduces inflammation, removes toxins, helps restore the integrity of the skin of the feet and nails, and improves their condition.</p> <p><b>Application:</b> daily for a month, in conjunction with the corresponding programs on DeVita Ritm (see "Clean skin" complex).</p>	No onychomycosis Without trichophyton Nails without fungus Without microsporium Without toxins	10:41

Healthy nails - automatic	It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month. In order to prevent the occurrence of the disease - once or twice a year (fall-spring).	No onychomycosis Without trichophyton Nails without fungus Without microsporium Without toxins	10:41
Stop skin fungus -automatic	The program inhibits the development of fungal flora on the scalp and body, reduces inflammation, reduces toxic load, and promotes a better quality of life. <b>Application:</b> daily for a month, in conjunction with the corresponding programs on DeVita Ritm (see "Clean skin" complex). It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month.	Without Epidermophyton Without trichophyton Without microsporium Fungus-free skin Without toxins	8:23
Low back without pain -automatic	The program relieves pain and improves blood flow in the lower back area. <b>Application:</b> daily for 2-4 weeks or as needed, in conjunction with the programs "Joints without pain-automatic", "Healthy spine-automatic" on DeVita Ritm .	Low back without inflammation Stop pain (low back)	6:14
Healthy back -automatic	The program helps to reduce muscle spasm, reduce pain, promotes muscle fiber repair, increases muscle elasticity, reduces muscle tone, and improves the condition of all tissues of the spine. <b>Application:</b> daily for a month, in conjunction with the programs "Joints without pain-automatic", "Healthy spine - automatic" on DeVita Ritm .	Straight spine Stop spasm (muscles) Muscles without tension	8:24



[www.deholding.info](http://www.deholding.info)

Deta-Elis Europa GmbH  
Kaiserstrasse 74, 63065 Offenbach, Germany