



Deta
Elis
Holding
Health
Innovation
Human

List of
manual programs

DeVita AP+
Model Mini

Description and application

N° p/p	The name of the program	Program description
Manual programs		
1	No staphylococcus	<p>According to statistics, 20% people are regular carriers of staphylococcus. Staphylococcus is a fairly common causative agent of infectious diseases affecting the skin, soft tissues, upper respiratory tract, and also the musculoskeletal system. The main cause of infection is often the hospital and doctors who become chronic carriers of bacteria and people who have taken antibiotics. Staphylococcus cause secondary diseases with influenza, wound infections, postoperative suppuration, and can also cause food poisoning. The program removes the cause of the infection, which can have a detrimental effect on the human body. It is also effective as a complex treatment for diseases caused by these bacteria.</p> <p>Application: daily up to 2 weeks, depending on the reduction of infection.</p>
2	No streptococcus	<p>Streptococcus is a widespread group of microorganisms. They usually live in the skin of a person, in the respiratory tract, intestines. They can cause acute inflammation of the ENT organs, urinary system, organs of the abdominal cavity, skin and soft tissues. Streptococcal infection is more severe than staphylococcal infection, which can lead to the appearance of joint, heart, kidney, sinusitis, etc. The program is effective for complex treatment, and separately.</p> <p>Application: daily up to 2 weeks, depending on the reduction of infection.</p>
3	No chlamydia	<p>Chlamydia cause sexually transmitted infections. Chlamydia often very asymptomatic, or with minimal manifestations. The organs of the genitourinary sphere are mainly affected. Women are more susceptible to this infection, children can become infected when passing through the birth canal of the mother. Chlamydia can cause serious complications, therefore it is very important to take timely measures to prevent the disease, observing hygiene rules and using DeVita technology</p> <p>Application: daily for 2-3 months.</p>

4	No helicobacter	<p>Helicobacter Piloni is a microorganism that lives in the human stomach and duodenum and provokes the development of many dangerous diseases (gastritis, ulcers, hepatitis and many others). According to the WHO, helicobacteriosis infected 60 of the world's population. Infection occurs very easily and is often considered a family disease, because if one member of a family becomes infected, the probability of infection of the rest is 95. Often the transmission of the parasite occurs as a result of the general use of personal hygiene products, utensils, with a kiss. Symptoms such as heartburn, bloating, pain, and heaviness in the stomach can be caused by this microorganism.</p> <p>Application: daily for 2 weeks.</p>
5	Stop bacteria	<p>Bacterial infections include a huge group of diseases caused by microorganisms - bacteria. The most important in the development of diseases caused by a bacterial infection is the release by the bacteria of toxins (poisons), which have an inflammatory effect, cause intoxication of the body and damage internal organs. The Stop bacteria program is used to reduce the severity of intoxication symptoms and alleviate the human condition during bacterial infections of the digestive tract and respiratory system caused by such forms as pneumococcus, listeria, and Klebsiella.</p> <p>Application: daily until the disappearance of signs of infection.</p>
6	No hepatitis A	<p>In developing countries and in places where hygiene rules are neglected, the incidence of viral hepatitis A is very high. Since the virus is transmitted through dirty hands, food, water, even salty. If you plan to travel to Asia and Africa, you need to take appropriate safety measures: observe hygiene, drink only bottled water, use the "Hepatitis A-free" program at the first suspicion of the presence of infection in you or others.</p> <p>Application: daily up to 2-3 times a day, at least 1 month.</p>
7	Without toxins	<p>The program promotes the removal of toxins, improves drainage functions, as a result, immunity and vitality are increased. Used as a detox program after parasitic, as well as independently for a month.</p>

7	Without toxins	<p>Unlike other detoxification programs, it does not expend energy resources of the body.</p> <p>Application: once a day after antiparasitic programs. Included in most automated programs.</p>
8	No herpes zoster	<p>The disease caused by this virus is called herpes zoster. If a person has had chickenpox (one of the forms of the zoster virus), the virus remains in the body in a latent form and hides in the nerve cells. With hypothermia, colds or chronic diseases, the virus comes to the surface of the skin with the formation of characteristic rashes on it and strong pain along any sensory nerve. The program aims to destroy the virus and remove the main signs of its manifestation.</p> <p>Application: daily up to 3-4 times a day until the disappearance of undesirable manifestations. At least 1 month. In conjunction with the "No herpes" program.</p> <p>Attention: the program is not used for chickenpox.</p>
9	No herpes simplex	<p>Herpes simplex has been known to mankind since ancient times. Every third person can be a virus carrier. It is localized on the mucous membranes of the lips, mouth, genitals. A decrease in immunity during hypothermia or overheating, during infectious diseases or during the period of menstruation in women contributes to the appearance of unpleasant rashes on the mucous membranes. The virus is transmitted primarily through contact with body fluids of an infected person (blood, saliva, sperm, mucous secretions). The program effectively reduces herpetic burden and prevents the appearance of signs of infection on the mucous membranes.</p> <p>Application: at the first sign of a rash every 3 hours, and then 1-2 times a day for 1 month. In conjunction with the program "No herpes"..</p>
10	No herpes	<p>Virus carriers are more than 90% of all people on the planet. You can get infected in various ways: through saliva, sex, hygiene items, as well as by airborne droplets. At the beginning, the disease can be asymptomatic, then characteristic signs appear: headache, malaise, fever, herpetic eruptions. The program will help reduce recovery time and reduce inflammatory and pain manifestations.</p>

10	No herpes	Application: 1-2 times a day for a month in conjunction with the No Herpes Zoster or No Herpes Simplex program.
11	Healthy lymph	<p>Lymph is the liquid tissue of the body. It plays a very important role in maintaining water balance, homeostasis of the microenvironment of cells and capable of leaching pathological microorganisms and toxins from tissues. In case of insufficiency of lymphatic circulation, the lymph does not cope with its main tasks, which leads to the emergence of a pathological process. The program helps to normalize lymphatic outflow, strengthen the walls of lymphatic vessels, improve the condition of lymph nodes, which leads to a decrease in swelling and better discharge of metabolites from the tissues. When using the program, the body's defenses increase.</p> <p>Application: if necessary, use 1 time per day or according to the schemes of the complexes (included in the "Cleansing" auto program).</p>
12	Anti-Cold	<p>The common cold is one of the most common infectious diseases in the world. The causative agents are respiratory viruses (adenoviruses, respiratory syncytial virus) that affect the upper respiratory tract (nose, pharynx, trachea, larynx). Diseases caused by this group of viruses are contagious, because they are transmitted by airborne droplets and are not treated with antibiotics. Complications after a cold - rhinitis, sinusitis, otitis media and more. Very often, an attenuation of the weakened organism can result in the accumulation of a bacterial infection. Therefore, it is advisable to use this program in combination with the antibacterial "No Staphylococcus", "No Streptococcus". The program is aimed at reducing inflammation of the mucous membranes of the nose, throat, bronchi and activation of local immunity.</p> <p>Application: at the first sign of a cold the first 3 days - up to 2-3 times a day. Then 1 time per day for 7-14 days</p>
13	No mould & fungi	<p>Microscopic fungi are ubiquitous. A modern heated apartment or private house is an excellent place for mold and fungus. Favorite mold habitats are bathrooms, libraries, and poorly ventilated living spaces.</p>

13	No mould & fungi	<p>Mushrooms can enter the body with food, water, air, with old dusty things. Mold fungi cause allergies, not only inhalation, but also a cross-reaction on the skin. Mold can penetrate the skin barrier wherever we bathe, where we walk barefoot, wherever there are conditions for mold (humidity, temperature, organic residues in the form of dirt or skin fat). Fungi and mold cause diseases of the respiratory organs, gastrointestinal tract, skin with a change in the tissue of these organs.</p> <p>To prevent fungal infections and reduce the occurrence of allergic reactions, it is recommended to use the Stop Mold and Fungi program. The program is effective against aspergillus.</p> <p>Application: 1-3 times a day for 4 weeks. For problems with skin and nails, they are used in a complex with the "Without fungus" program.</p>
14	No candida	<p>The program protects you from infection with such a common fungus as candida. It is part of the normal microflora of the intestine, vagina, and mouth. But with a decrease in immunity, candida begins to multiply actively and harm the body. Diseases are also caused by pathogenic candida, which enter the body with food, through household items, after the use of potent antibiotics, during intercourse with an infected partner. Symptoms are manifested in the form of allergic rashes on the skin, thrush, pain and burning in the genitals, dysbiosis. Prevent candida breeding using the "No candida" program.</p> <p>Application: the program should be applied 1-3 times a day for 2-4 weeks. Complete with the program "No to the simplest".</p>
15	No fungi	<p>Human fungal diseases belong to the group of infectious diseases caused by pathogenic fungi. The main danger of fungal diseases is that often the areas affected by the fungus are an open road for pathogens of many different infections. Most often, mycoses affect the skin, its derivatives (hair, nails), mucous membranes of internal organs.</p> <p>Treatment requires an integrated approach and hygiene. The program actively fights against such fungi as trichophyton, microsporium, and mucor.</p>

15	No fungi	<p>Application: In order to prevent infection with fungal infections, it is recommended to use the program every time after visiting the sauna, pool, beach, communicating with an infected person. If there is a fungal infection in the body, the program is used 1-2 times a day for a month.</p>
16	No protists	<p>The program protects you from the negative effects of a large number of protozoa. Every fifth person suffers from them, and he becomes infected through dirty hands or affected food and water. The infected person experiences severe abdominal pain, bloating, constant feeling of fatigue, depression. The liver and pancreas may be affected. Of course, all these signs interfere with your life, and besides, you can infect your loved ones. The program works on the amoeba, lamblia, destroying them. Observing the rules of personal hygiene, as well as the "golden rule" - not to drink tap water and regularly using the program "No to the simplest" as a preventive measure, you will be protected from these microorganisms.</p> <p>Application: it is recommended to use the program 1-2 times a day for 10-14 days. With re-infection, it is necessary to identify the sources of infection.</p>
17	Intimate health	<p>The causative agents of diseases of the genitourinary system are many and often the cause of the disease is not one microorganism, but several. Infections can occur in acute or chronic form. Chronic infection over time can cause narrowing of the urethra, inflammation of the joints, lead to infertility and increase the risk of an ectopic pregnancy. Periodically, during a cold or emotional overload, exacerbations may occur. Infections are mainly transmitted through sexual intercourse, and some by birth from a sick mother. Parasites such as ureaplasma can live in the body for years without causing any symptoms and not manifesting themselves. The main factor in protecting the body is the normal microflora, and modern DeVita technologies that adversely affect the parasitic flora, while maintaining the body's natural defenses. The program effectively fights against parasites such as ureaplasma, toxoplasma, mycoplasma.</p> <p>Application: daily for a month.</p>

18	No schistosomes (blood-flukes) & trematods	<p>Trematodes are flatworms that can cause a variety of serious diseases in humans. They can parasitize in the blood (schistosomes), liver, intestines, lungs. Infection occurs when chewing blades of grass, especially near ponds or while bathing in them, when eating poorly washed greens or lettuce, and also not fully cooked cattle meat, fish. The program is aimed at eliminating most types of trematodes, including schistosomes, regardless of the stage of their development and localization in the body.</p> <p>Application: daily for 14 days.</p>
19	No nematodes and ascarides	<p>With unwashed fruits, poorly prepared meat, and dirty hands, worm eggs that cannot be seen with the naked eye can get into our mouths. Worms, living in the host's body, feed on its tissues, "take food" from a person (that's why you want to eat a lot), and release toxic substances. A person is losing weight, his hair is falling out, his skin is peeling. The most common helminths are roundworms. They can cause damage to the wall of the small intestine, intestinal obstruction, inflammation of the peritoneum, bronchitis, focal pneumonia, hemoptysis, liver and pancreas diseases. Ascaris parasitism causes an allergic reaction, which can manifest itself in the form of bronchial asthma and skin rash. Under the action of the program, worms cease to feed, which leads to their death and the natural removal of parasites from the digestive tract. The program does not have a toxic effect on the human body. It is recommended that people with pets, small children and occupations related to working with animals regularly conduct preventive measures once a quarter.</p> <p>Application: the program is used daily. A total of at least 10 sessions. After 14 days, repeat the course.</p>
20	No parasites, General	<p>Parasites are infesting the human body through the skin, and can also be carried into the body with food, water, and dust. Many people who lead a healthy lifestyle experience health problems due to parasites. Improving the body (proper nutrition, exercise, hardening procedures) without ridding the body of parasites does not produce a pronounced positive effect. It is parasites that are the trigger for many chronic diseases: cholecystitis, cholelithiasis, pancreatitis, colitis, diabetes mellitus, bronchial asthma, and dermatitis.</p>

20	No parasites, General	<p>Problems such as chronic fatigue, irritability and anxiety, hyperactivity in children, anemia, brittle nails and hair, problem skin, headaches, appetite disorders, decreased immunity can be signals of ongoing parasitosis. All parasitic diseases are terrible in that they perfectly mask for other diseases. Prevention of parasitic diseases is, first of all, observance of hygienic norms of social behavior, use of the "Stop parasites common" programs. These measures are aimed at maintaining your health and the health of others.</p> <p>Application: the program is used daily. A total of at least 10 sessions. After 14 days, repeat the course.</p>
21	No parasites, Digestive tract	<p>The presence of parasites in the digestive tract and liver causes the development of diseases such as cholangitis, biliary dyskinesia, cholecystitis, hepatitis. Their presence in the intestine causes inflammation of the intestinal walls and a decrease in the absorption of nutrients. Toxic substances that secrete parasites cause pain in the joints and muscles, affect the central nervous system. Contribute to the onset of depressive states, nervousness, anxiety, and sleep disturbances. The program is aimed at cleansing the body from parasites of the gastrointestinal tract, liver and pancreas.</p> <p>Application: the program is used daily for 7-10 days. After 14 days, repeat the course.</p>
22	Anti-Flu	<p>Influenza is a severe viral infection that affects people of any age. Up to 500 million people get the flu every year. In severe influenza, irreversible damage to the cardiovascular system, respiratory tract, and central nervous system occurs. The virus is transmitted very easily by airborne droplets (by sneezing, coughing, talking). The most common complication of influenza is pneumonia. The program is aimed at the destruction of various strains of influenza viruses, the removal of inflammatory manifestations.</p> <p>Application: to prevent influenza during an epidemic, the program is applied once every 2 days. If necessary - daily up to 2-3 times a day, alternating it with the Stop Flu Basic program. After that it is necessary to use the "General drainage" program. It can be used together with the programs "No to staphylococci", "No to streptococci", "No candida" in order to avoid bacterial and fungal burden, as well as the programs "Anti-virus", "Anti-cold".</p>

23	Against flu - General	<p>Influenza is a disease that affects the upper and lower respiratory tract, is accompanied by severe intoxication and can lead to serious complications. Influenza always begins acutely and quickly (48 hours pass from the moment of infection to the first signs). The main symptom of intoxication is an increased body temperature (40 °C and higher). This temperature can persist for more than 5 days, which causes: shortness of breath or respiratory failure, impaired consciousness, convulsions. In people with poor health, elderly, young children, this temperature leads to an infectious toxic shock (a condition when vital organs fail). Using the program "Stop influenza basic" you can alleviate the patient's condition, and with the integrated use of several programs against influenza viruses, you can quickly restore your health without complications.</p> <p>Application: at the first sign of flu with increased body temperature, use the program 1-3 times in a row or alternating it with the Anti-flu program. Be sure to drink plenty of alcohol and observe the rules of personal hygiene.</p>
24	Healthy gastrointestinal	<p>Intestinal infection - a disease that occurs when certain pathogens of viruses or bacteria enter the body. People of all ages are affected by this infection, most often children and those who take antibiotics. An intestinal infection is accompanied: nausea, vomiting, diarrhea, abdominal pain, fever. In case of viral infections, other symptoms can also be observed: runny nose and redness of the throat, as well as pain that occurs when swallowing, lack of appetite, loss of strength. The program is aimed at the destruction of the intestinal flu virus - rotavirus and pathogenic strains of E. coli, will help normalize the condition and stop dehydration destructive for the body.</p> <p>Application: Use on the first day 4-5 times in a row to achieve a stable result. Then within 3-5 days, 1-2 times daily, alternate with the program "General drainage" and "Anti-virus". Be sure to drink as much water as possible.</p>
25	Help with poisoning	<p>When ingested with food, water, enterobacteria or enteroviruses through the mucous membranes of the respiratory tract, an infection develops, which is called poisoning in the people. This infection can have a latent character or obvious signs of its presence: headache, muscle pain, in the abdominal region, fever, vomiting, upset stool. Subsequently, it can lead to serious organ damage.</p>

25	Help with poisoning	<p>The program helps the body cope with the infection caused by enteroviruses and some types of salmonella, accelerates the recovery process, and reduces the likelihood of complications.</p> <p>Application: Use on the first day 4-5 times in a row to achieve a sustainable result. Then, 3-5 days daily 1-2 times, alternate with the program "General drainage" and "Antivirus". Be sure to drink as much water as possible.</p>
26	Antivirus	<p>Viral infections are a large group of diseases caused by the virus. Viruses are a transitional form between animate and inanimate nature. Since viruses have genetic material, but do not have their own cell, they are only intracellular parasites. This is precisely the difficulty of their destruction. Viruses can be transmitted by airborne droplets, through various objects, kisses and sexually. A large role in the complete extermination of viruses in the body belongs to the immune system. The Anti-Virus program is aimed at combating various viruses, and in particular the cytomegalovirus and Epstein-Barr virus. It is also used to attenuate any strains of viruses and enhance the body's immunity.</p> <p>Application: daily for a month.</p>
27	Drainage - General	<p>One of the main programs to cleanse the body of toxins. The lymphatic system is involved in the removal of toxic substances from tissues and organs. With various diseases, the amount of harmful substances increases in many times. The lymphatic system does not have enough time to remove them and intoxication of the body may occur again. In addition, due to illness or injury, a violation of the natural outflow of contents from the wound cavity or organ is possible. Excessive accumulation of fluid containing pathogens can cause serious inflammation and deterioration of the human body. Using the "Drainage - General" program you protect your body from the excessive toxic load that it experiences with various diseases or injuries. The program is aimed at neutralizing and removing toxins from the body, while not loading the lymphatic system.</p> <p>Application: once after using any of the DeVita AP device programs. Included in the "Cleanse" auto-program.</p>

28	Drainage Antitox	<p>The intercellular matrix is the main place of accumulation of toxins, which leads to metabolic disorders, the development of many serious chronic diseases, and premature aging of the body. Tissue drainage, is a means of detoxification of the tissue sector, helps to remove salts of heavy metals and radionuclides, improve metabolic processes and the general condition of the body.</p> <p>Application: with severe intoxication daily and in conjunction with the "General drainage" program.</p>
29	Mexican detoxification	<p>The system of deep complete cleansing of the body, which was proposed by R. Rife. The program helps to eliminate metabolic products and toxins, improve the drainage function of the lymphatic system, increase the detoxification capabilities of the liver, activate the excretory function of the kidneys, and increase the level of immune defense.</p> <p>Application: daily for a month in the absence of serious diseases of the liver, kidneys, lymphatic system. If there are any, it is recommended to use the automatic program "Complete Cleansing" DeVita Ritm for gentle cleansing of the body.</p>
30	No allergies	<p>Allergy is an increased sensitivity of the body to any substances - products, wool, dust, pollen, chemical, infectious agents. Currently, more than 85% of the entire population of our planet suffer from its various varieties. The cause of allergies is most often malnutrition and an unhealthy lifestyle. An allergy can be recognized by an unexpected runny nose, sneezing or lacrimation, redness and itching of the skin. The program removes the main manifestations of allergies, eases the human condition, and helps stabilize the immune system.</p> <p>Application: daily if indicated and to alleviate the condition.</p>
31	Weight Control	<p>The program aims to combat adenovirus 36, a common cause of overweight. In addition, when it is used, the ANS is regulated, and the processes of weight loss are launched. The program also helps maintain the body's energy resources.</p> <p>Application: it is advisable to use according to the scheme of the complex "Slim figure".</p>

32	Without Salmonella	<p>Salmonella is a representative of the enterobacteriaceae family. Resistant to the external environment: can be stored in house dust, in open water, in frozen meat and sausages, poultry eggs, dairy products from 3 to 6 months. They perish at a temperature of more than 70 degrees. Infection is accompanied by: high fever, abdominal pain, vomiting, loose stools. However, carriage may be without symptoms.</p> <p>Application: daily for 2 weeks.</p>
33	Without Yersinia	<p>One of the representatives of enterobacteria, which causes inflammation in the small intestine, but can also affect other organs: liver, kidneys, heart, joints. It is found in water, soil, on food. It tolerates low temperatures, so it can be stored and multiplied in the refrigerator. In the human body, without causing signs of inflammation, it can be located until the immune system is weakened. Then the exacerbation begins.</p> <p>Application: daily for 2 weeks.</p>
34	Without Proteus	<p>Bacteria of the genus Proteus are part of the conditionally pathogenic intestinal flora. With active growth, they cause damage to the digestive tract, urinary tract, congenital malformations, and wound infections. The pathogenicity of bacteria is observed in the case of uncontrolled intake of antibiotics. Or when bacteria enter from the outside through food, water, unwashed hands.</p> <p>Application: daily for 2 weeks.</p>
35	Without Klebsiella	<p>Conditionally pathogenic bacterium. With an increase in the number, it affects the organs of the gastrointestinal tract, respiratory tract, and urinary system. Causes severe intoxication with impaired functionality of these organs. Ways of infection through food, water, when communicating with an infected person.</p> <p>Application: daily for 2 weeks.</p>
36	Without Campylobacter	<p>Campylobacter bacteria are causative agents of intestinal infection. Introduce into the mucosa and lymphatic tract of the intestine and other organs. Infection occurs through agricultural animals and birds, through the use of contaminated meat products, milk, water. Less often through damaged skin, with bites of infected animals. Also transplacental or during childbirth. At risk are rural residents, livestock and poultry farmers, as well as tourists visiting developing countries.</p> <p>Application: daily for 2 weeks.</p>

37	Without Clostridia	<p>Clostridium difficile is part of the normal intestinal microflora, but with dysbiosis can cause damage to the walls of the colon. Clostridium botulinum is an inhabitant of soil and sea water, distributed everywhere and form the most powerful bacterial toxin known - botulinum toxin.</p> <p>Application: daily for 2 weeks.</p>
38	Without Shigella	<p>Bacteria that cause dysentery. Infection can occur by fecal-oral or contact-household route, through water, food. Carriers can be flies and cockroaches.</p> <p>Application: daily for 2 weeks.</p>
39	Without Rotavirus	<p>The intestinal flu virus multiplies in the cells of the gastrointestinal mucosa, which leads to severe diarrhea and dehydration of the body. The main route of transmission is food.</p> <p>Application: daily for 2 weeks.</p>
40	Anti-virus Gastrointestinal	<p>Enterovirus and Koksaki virus belong to the intestinal group, are able to affect the central nervous system, muscles, heart, skin, respiratory and genitourinary systems. They are transmitted from person to person by airborne droplets and in everyday life. The wave of infection is in the warm season.</p> <p>Application: daily for 2 weeks.</p>
41	Without Trematodes - General 1	<p>Trematodes, or flukes, are parasitic flatworms that infect the liver, pancreas, intestines, lungs, blood vessels, and brain.</p>
42	Without Trematodes - General 2	<p>Ways of infection: swimming in open waters of countries of Asia, Africa, South America, eating poorly cooked or raw fish, shellfish, crustaceans, and meat.</p> <p>Application: programs are used together daily for 2 weeks.</p>
43	Stop Amoeba	<p>Entamoeba histolytica is a representative of protozoa. It can cause amoebic dysentery, accompanied by damage to the colon, liver, lungs and other organs. It is ingested with food and water.</p> <p>Application: daily for 2 weeks.</p>

44	Without Enterococcus	<p>Enterococcus – present normal human microflora. When using antibiotics, impaired barrier function of the skin and mucous membranes, serious diseases and in older people, enterococcus can cause problems in the small intestine, urinary tract, vagina, heart, liver, skin.</p> <p>Application: daily for 2 weeks.</p>
45	Without nematodes and ascarididae - New	<p>A new program to combat round parasitic worms. Covers the majority of nematode species that infect humans and animals.</p> <p>Application: daily for 2 weeks.</p>
46	Stop Adenovirus - General	<p>Adenoviruses are a group of viruses that cause acute diseases, with damage to the respiratory system, eyes, and lymph nodes. Infection occurs by airborne droplets, often in the cold season.</p> <p>Application: daily for 2 weeks.</p>
47	Stop Coronavirus	<p>The virus is spread in Europe, Asia, North America. It is transmitted by airborne droplets with damage to the respiratory tract, up to the destruction of the pulmonary alveoli.</p> <p>Application: daily for 2 weeks.</p>
48	Stop Penicillium	<p>Genus of molds affecting the organs of the urinary and respiratory systems, lymph nodes, skin. They cause severe intoxication, are allergens.</p> <p>Application: daily for 2 weeks.</p>
49	Without Rickettsia	<p>Intracellular parasite. It is able to affect any human organs, cause significant disorders of the central nervous system and circulatory disorders. Routes of transmission: through a tick bite, food, contacts. Distributed in developing countries.</p> <p>Application: daily for 2 weeks.</p>
50	Without Rickettsia	<p>HPV affects the skin and mucous membranes with the formation of growths. Some can cause malignant tumors, in the risk zone of women. Ways of transmission of this infection: contact-household, sexual, in childbirth from mother to child.</p> <p>Application: daily for 2 weeks.</p>

51	Without Papilloma	<p>With the defeat of HPV of the mucous surface of the urinary system, there are formations in the renal pelvis and bladder. With prolonged exposure to adverse factors, these formations can not only reduce the functionality of organs, but also are capable of degeneration.</p> <p>Application: daily for 2 weeks.</p>
52	Without Trichomonas	<p>Trichomonas is one of the simplest pathogens of many sexual diseases. Infection occurs with any type of sexual contact, as well as through household contact: the rim of the toilet bowl, the pool, the use of common towels, washcloths, hygiene products, and fitting other people's clothes (underwear). Once in the bloodstream, they can disguise themselves as blood cells. The biggest threat, most often, is not in the simplest, but in the infectious agents he carries: chlamydia, candida, gonococcus, cytomegalovirus and herpes viruses.</p> <p>Application: daily for 2 weeks.</p>
53	Without Leptospira	<p>A bacterium that causes damage to capillaries, liver, kidneys, and muscles. Distributed in all regions. Infection occurs from domestic and wild animals, through damaged mucous membranes and skin, from natural water sources.</p> <p>Application: daily for 2 weeks.</p>
54	Without Cryptococcus	<p>Yeast fungi of the genus Cryptococcus, affect the central nervous system, lungs, skin, mucous membranes. The causative agent is present in soil contaminated with bird droppings. The transmission path is dusty.</p> <p>Application: daily for 2 weeks.</p>
55	Without Toxoplasma	<p>The prevalence of Toxoplasma infection in the world is very high up to 90. Being in the body often does not cause problems for its owner. Infection occurs either in utero or is acquired. Toxoplasma affects the nervous system, eyes, heart, muscles, lungs, intestines, skin.</p> <p>Application: daily for 2 weeks.</p>
56	Without Ureaplasma	<p>Ureaplasma is a pathogenic form that causes diseases of the genitourinary system, joints, pathology of pregnancy and the fetus. The transmission of an infectious agent occurs through sexual contact, in utero and during childbirth.</p> <p>Application: daily for 2 weeks.</p>

57	Without Echinococcus	<p>A serious problem is represented by tapeworm larvae, which form cysts in the body. These formations compress the surrounding tissue, which leads to atrophy and the death of the latter. Larvae affect the liver, lungs, kidneys, bones, and brain. Infection occurs through contact with invasive animals, when collecting berries and herbs, drinking water from infected sources. Also at risk are slaughterhouse workers, shepherds, tanners.</p> <p>Application: daily for 2 weeks.</p>
58	Without Pseudomonas aeruginosa	<p>Pseudomonas aeruginosa is a bacterium characterized by natural resistance to most antimicrobial agents. It affects the organs of the gastrointestinal tract, respiratory system, central nervous system, other organs and systems. The source is people and animals, patients and carriers. Transmission ways: contact-household, airborne, food.</p> <p>Application: daily for 2 weeks.</p>
59	Against Lyme Disease 1	<p>Tick-borne borreliosis (Lyme disease) is the most common problem that occurs when a person bites an infected ixodid tick. The central nervous system, heart, joints are affected.</p> <p>Application: the program is used together with Stop Ticks 2, after bites every day for 2 weeks. For problems arising from bites in the past, see the SOS Pliers complex.</p>
60	Against Lyme Disease 2	<p>Tick-borne encephalitis virus is transmitted by bites of ixodid ticks, less often through the use of raw milk from cows and goats infected with the virus. Distribution: Siberia, the Far East, the Urals. The main organ of the central nervous system.</p> <p>Application: the program is used together with Stop Ticks 1, after bites every day for 2 weeks. For problems arising from bites in the past, see the SOS Pliers complex.</p>
61	Healthy Cell - 1	<p>Various pathogenic factors contribute to the accumulation of toxins in the body that can cause cell degeneration. The programs "Without mutations 1 and 2" reduce the risk of oncological pathology, contribute to the correction of pre-tumor changes, reduce inflammatory processes and enhance immune defense.</p> <p>Application: programs are used together, see the "Protect Yourself" complex.</p>
62	Healthy Cell - 2	

63	Healthy Gums	<p>Gum disease (gingivitis) is caused by viruses, bacteria, and fungi. At risk may be people wearing dentures and braces. The program helps to reduce the pathogenic load, reduce the inflammatory process, and improve trophic tissue.</p> <p>Application: daily for 1-2 months.</p>
64	Mucous Membrane Regeneration	<p>The appearance of ulcers (aphthae) on the oral mucosa indicates the presence of an inflammatory process that can be caused by viral or bacterial agents, trauma or a lack of trace elements. The program removes the inflammatory process, helps restore the oral mucosa, and strengthens local immunity.</p> <p>Application: daily for 2 weeks. It is advisable to use it in conjunction with the "Mouth in order" program.</p>
65	Healthy Mouth	<p>Inflammation of the oral mucosa (stomatitis) occurs due to a number of factors: the presence of infectious agents, damage to the oral cavity, during dental interventions, and an unbalanced diet. It is localized on the tongue, the inner surface of the lips, cheeks, and palate. The program removes the inflammatory process, helps restore the oral mucosa, and strengthens local immunity.</p> <p>Application: daily for 2 weeks. Can be combined with the program "Renewal of the mucosa". Use in any dental procedures accompanied by mucosal damage.</p>
66	Without Fistula	<p>Suppuration of the gums is very often accompanied by the formation of a fistula. This can be fraught with serious consequences from tooth loss to deep tissue inflammation. The program helps prevent suppuration, reduces the inflammatory process.</p> <p>Application: daily for 2-4 weeks. Use after dental treatment or tooth extraction.</p>
67	Stop Injury (teeth)	<p>The program reduces pain, speeds up the healing process, and relieves swelling.</p> <p>Application: daily for 1st-2nd week. Use immediately after dental procedures for the treatment, removal, implantation of teeth.</p>
68	CNS Without Problems	<p>The program helps to improve the functioning of the nervous system after exposure to various factors (infectious, traumatic, etc.).</p>

68	CNS Without Problems	<p>It helps to restore nerve connections between the central nervous system and peripheral, relieves swelling and pain.</p> <p>Application: for 1-6 months daily. Preferably together with the program "Nerves are in order" (see the SOS Ticks complex).</p>
69	Healthy Nerves	<p>In case of violation of the integrity of the nerve fiber, functional disorders of the nervous system can occur. The program prevents the destruction of the myelin sheath of nerve fibers, thereby preserving the conductivity of the nerve fiber.</p> <p>Application: for 1-6 months daily. It is advisable with the program "central nervous system without problems" (see the complex "SOS Ticks").</p>
70	Clean Sinuses	<p>Inflammation of the mucous membrane of one or more sinuses - sinusitis. It can occur as a complication of infectious diseases of the nose and respiratory tract, as well as after facial trauma.</p> <p>The program helps to reduce the inflammatory process in the sinus area, better discharge of mucus, and reduce pain during acute processes.</p> <p>Application: daily for 2 weeks.</p>
71	Clean Sinuses - General	<p>The greatest problem for many people is chronic sinusitis, which occurs without obvious manifestations, but for a sufficiently long time. The causes may be infections of the upper respiratory tract, the presence of polyps, curvature of the nasal septum, and carious teeth. The program helps to reduce the inflammatory process in the sinus area, reduce swelling, and prevent exacerbations.</p> <p>Application: daily for a month if there is a problem. And in order to prevent exacerbation during the cold season.</p>
72	Ears Without Inflammation	<p>The causes of ear inflammation can be: a pathological process in the area of the nasopharynx, curvature of the nasal septum, ear injury.</p> <p>The program helps to reduce swelling, reduce inflammation and pain.</p> <p>Application: daily for 2 weeks. Alternate with DeVita Ritm Healthy Ears.</p>

73	Throat Without Pain	<p>Tonsil inflammation of the pharyngeal ring is often observed after suffering tonsillitis or other infectious diseases of the pharynx. May occur in acute or chronic form. Transmission externally - by airborne droplets, or internally - from internal foci of chronic infection. The program helps reduce tonsil inflammation, swelling, pain, and toxic stress.</p> <p>Application: daily for 2 weeks. Alternate with DeVita Ritm "Throat no problem".</p>
74	Against Cold - General	<p>Inflammation of the mucous membranes of the nasopharynx (nasopharyngitis) has an infectious nature, viral or bacterial. The highest frequency is observed in the offseason. The program reduces inflammation, swelling, prevents the development of exacerbations.</p> <p>Application: daily for 2 weeks.</p>
75	Anti-inflammation	<p>Inflammation is a body reaction that occurs in response to damage or the action of a pathogenic agent. The program helps the body eliminate the source of inflammation, wherever it is located, reduce the influence of pathogenic factors, and also restore homeostasis.</p> <p>Application: daily for 1-4 weeks. It is advisable to combine with any other programs in order to reduce the inflammatory process.</p>
76	Stop Migraine	<p>Migraine headaches are one of the most common types of headaches. It is characterized by acute throbbing pain localized in a certain part of the head. The program relieves spasms of blood vessels, reduces pain, helps normalize blood circulation in the head.</p> <p>Application: as needed.</p>
77	Antiseptic	<p>Damage to the skin: scratches, cuts, burns, etc. are open gates for infectious agents. As a result, the wounds may not heal for a long time, become inflamed, or hurt. The program helps to disinfect the wound, reduce inflammation and quickly tighten it.</p> <p>Application: daily for 1-4 weeks. It is advisable to combine with any other programs in order to reduce the inflammatory process.</p>

78	Healthy Bile	<p>Improper and irregular nutrition (prolonged fasting, alcohol, spices, overly spicy fatty and overcooked meals, snacks, refusal to consume vegetable oil), stress, a small amount of fluid drunk, infectious liver lesions are the causes of problems in the work of the gall the bubble. The program facilitates the outflow of bile, helps restore the work of the gallbladder, reduces pain.</p> <p>Application: daily for 1-6 months, preferably in the evening. Otherwise, in a day.</p>
79	Without cough	<p>Cough is a protective reflex reaction of the body that occurs due to irritation of the respiratory tract mucous membrane. It can be a sign of inflammation of the respiratory system, which can be caused by pathogenic viruses, bacteria, fungi and other agents, as well as allergens. The program alleviates the condition in the presence of acute or chronic cough, reduces inflammation of the mucous membrane of the respiratory tract.</p> <p>Application: daily for 2-4 weeks.</p>
80	Without bacterial pneumonia	<p>According to scientists, the first place among the causes of lower respiratory tract inflammation is occupied by microbial infections. Most bacteria that affect the respiratory system are part of the natural micro flora of the organism, but under adverse conditions they become pathogenic. The program fights against such microorganisms as: Streptococcus pneumoniae, Chlamydomphila pneumoniae, Mycoplasma pneumoniae, Klebsiella pneumoniae, etc.</p> <p>Application: daily for 2 weeks.</p>
81	No Coronavirus - New	<p>A new improved and expanded program that allows you to fight with all types of coronavirus, including new mutated forms.</p> <p>Application: daily for 2 weeks. In order to prevent the disease 2-3 times a week.</p>
82	Without pneumonia	<p>Inflammation of the lung tissue can be caused by pathogenic flora of various forms: viral, bacterial, fungal, etc. It arises not only as an independent form, but also as a result of the common cold, decreased immunity, and lack of vitamins. The program relieves from swelling and inflammation of the lower parts of the respiratory system, improves the functionality and energy of these organs, increases local immunity.</p> <p>Application: daily for 2-4 weeks. It is advised to combine with the program "Without cough".</p>

83	Low back without inflammation	<p>As a result of pinched and inflamed spinal nerve roots in the lumbar and/or sacral spine, many people may experience sharp pain, often spreading into the leg. The most susceptible to the disease are persons who, due to their occupation, are forced to lift heavy weights, be in a bending position or sit for long periods of time. The program relieves pain and reduces inflammation.</p> <p>Application: daily for 2-4weeks.</p>
84	Stop acne bacteria	<p>Propionibacterium is a conditionally pathogenic flora of human skin. Under the right conditions, as a result of blockage of the sebaceous gland ducts and increased sebum production, a favorable environment is created for their reproduction inside the hair follicles with the formation of acne (inflammation of the sebaceous glands). The activation of the pathogenetic properties of Propionibacteria may be linked to sex, age, and genetics. The program fights pathogenic microflora, and reduces skin inflammation.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" and "Skin without acne" programs on DeVita AP and the corresponding programs on DeVita Ritm (see Clear skin complex)</p> <p>In case of skin predisposition to acne formation, in order to prevent inflammation, it is desirable to use the program every quarter according to the above scheme.</p>
85	Without microsporium	<p>A fungus of the genus Microsporium primarily affects the skin, scalp, and hairy areas of the skin (beard, mustache), causing an infectious disease in animals and humans (ringworm). It is mainly transmitted from animals, less frequently from person to person or through household items. The program destroys pathogenic fungi, and it improves the condition of human skin.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see Clear skin complex).</p>
86	No onychomycosis	<p>Onychomycosis is a fungal lesion of the nails, usually caused by fungi of the genus Trichophyton, Epidermophyton floccosum and microsporidia.</p>

86	No onychomycosis	<p>Men and the elderly are most likely to be affected. The development of onychomycosis provoke trauma, flat feet, varicose veins, hyperhidrosis, endocrine and immunological disorders, blood diseases, vascular pathology, prolonged use of hormonal drugs, antibiotics and cytostatics, wearing tight shoes and refusal to comply with common norms of personal hygiene. The program causes the death of pathogenic fungi, reduces inflammation, and helps restore the integrity of the skin and nails.</p> <p>Application: daily for a month, in conjunction with the "Nails without fungus", "Without trichophyton", "Without toxins" programs on DeVita AP and the corresponding programs on DeVita Ritm (see "Clean skin" complex). It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month. In order to prevent the occurrence of the disease - once or twice a year (fall-spring).</p>
87	Zygomycetes, stop!	<p>People with low immune defenses are most likely to get infected with these fungi. Persons at risk include long-term users of immunosuppressive medications, persons with diabetes mellitus, cancer or acquired immunodeficiency syndrome. The fungus affects internal organs, more commonly the CNS, lungs and GI tract. The program inhibits the development of fungal flora, and it improves the person's condition.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm depending on the concerned organ. To prevent infection - 2 times a year for people with the above-mentioned health problems.</p>
88	Nails without fungus	<p>Mycosis of feet and nails is an infectious disease caused by pathogenic or conditionally pathogenic fungi of the species Trichophyton, Epidermophyton, Candida, Scytalidium. Risk factors are: violation of personal hygiene; sharing shoes; visiting public baths, swimming pools, and beaches; climatic conditions; injuries; wearing tight shoes; chronic diseases, etc.</p> <p>The program is aimed at destroying pathogenic flora, reducing inflammation, restoring the condition of the skin of the feet and nails.</p>

88	Nails without fungus	<p>Application: daily for a month, in conjunction with the "No onychomycosis", "Without trichophyton", "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see Clear skin complex). It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month. In order to prevent the occurrence of the disease - once or twice a year (fall-spring), in case of visiting a swimming pool, bath, beach - every time after the visit.</p>
89	Without trichophyton	<p>Fungi of the genus Trichophyton affect the scalp, smooth skin and nails. Infection occurs through contact with a sick person or their belongings, as well as through contact with animals and their care items. The disease can occur in both acute and chronic forms. The program inhibits the development of fungi and inflammation.</p> <p>Application: daily for a month, in conjunction with the "Without trichophyton", "Nails without fungus", "Without toxins" programs on DeVita AP and the corresponding programs on DeVita Ritm (see "Clean skin" complex). It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month. In order to prevent the occurrence of the disease - once or twice a year (fall-spring).</p>
90	Mycosis fungoides	<p>Mycosis fungoides affects primarily the lymphoid-reticular system of the skin, and subsequently the internal organs. The program helps to reduce the spread of fungus into the tissues and reduce the development of inflammation.</p> <p>Application: daily, in conjunction with the "Without toxins" program on the DeVita AP and the "Healthy day" and "Health" programs on the DeVita Ritm.</p>
91	Solar keratosis	<p>The main cause of formation is due to prolonged exposure of the skin to direct sun rays, but manifests itself only with age when the immune forces are weakened. Therefore, it is most common in the elderly. Factors such as heredity, wearing tight clothing, taking certain medications, etc. can be triggers of the disease. The program helps to reduce the chances of keratosis manifestation and progression.</p> <p>Application: daily, in conjunction with DeVita Ritm programs (see "Clean skin" complex).</p>

92	Stop warts	<p>Warts are benign skin neoplasms caused by human papillomaviruses. The transmission of papillomavirus can occur through household and sexual contact. The program fights HPV, preventing them from spreading in the body and forming neoplasms. It is necessary to consult a dermatologist regarding the removal of formations from the skin.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the programs DeVita Ritm (see Clear skin complex). Family members and sexual partners should use the program at the same time to prevent re-infection.</p>
93	Stop pimples	<p>Inflammation of hair follicles, provoked by the spread of staphylo-streptococcal infection with poor personal hygiene and/or the presence of chronic intoxication, infection, hypovitaminosis, leads to the formation of furuncles. The program inhibits the activity of bacterial flora, reduces inflammatory manifestations.</p> <p>Application: daily for a month, in conjunction with the program "Mode 2 - automatic".</p>
94	Fingers without inflammation	<p>In the presence of microtraumas and injuries to the perionychium area of the fingers or toes and the development of bacterial, most often - staphylococcal, less often - streptococcal or mixed flora in the tissues of the fingers can lead to acute purulent inflammation - paronychia. The program reduces tissue inflammation, swelling, improves overall health, and reduces the risk of complications.</p> <p>Application: daily for a month, in conjunction with the program "Mode 2 - automatic".</p>
95	No scleroderma	<p>Scleroderma is a connective tissue disease associated with impaired blood supply and subsequent thickening of tissues. The program improves blood supply to the tissues, slowing down the thickening processes.</p> <p>Application: daily for a month, in conjunction with the program DeVita Ritm. The cycle is repeated 2 times a year.</p>

96	Erysipelas, stop	<p>Erysipelas is an infectious disease affecting the skin and mucous membranes. It is transmitted by airborne droplets, sometimes by contact. Injuries and microtraumas of the skin and mucous membranes serve as entry gates for infection. The program prevents the spread of the disease, helps to reduce inflammation and intoxication.</p> <p>Application: daily, in conjunction with the "Without toxins", "No streptococcus" programs on the DeVita AP and the "Healthy day" and "Health" programs on the DeVita Ritm. The programs should be used in cycles 2-3 times a year for a month.</p>
97	Relieve psoriasis	<p>Psoriasis is a chronic non-infectious skin disease. Occurrence is associated with hereditary predisposition, immunity disorders, endocrine and nervous system disorders. The program helps reduce inflammation, prevent complications, and improve the quality of life.</p> <p>Application: daily, in conjunction with DeVita Ritm programs (see "Clean skin" complex).</p>
98	Leprosy	<p>Leprosy or leprosy is a chronic infectious disease affecting mainly the skin, peripheral nerves, upper respiratory mucosa and visual organs. It is characterized by a prolonged cycle (from a year to decades) with periodic exacerbations and can be disguised as other diseases, most often skin diseases. The program is aimed at destroying bacteria - the cause of the disease, reducing the inflammatory process, and preventing complications.</p> <p>Application: in conjunction with the "Without toxins", "Restoration 2- Automatic," programs on the DeVita AP and the "Healthy day" and "Health" programs on the DeVita Ritm.</p>
99	Without Epidermophyton	<p>A genus of fungi that cause lesions (mycoses) of the skin, hair, and nails. It is the most common group of mycoses and occurs in people of all ages, sexes, occupations and health conditions. More commonly settles in the folds of the skin. The program inhibits the development of fungal flora, reduces inflammation, and improves the quality of life.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see Clear skin complex). It is mandatory to maintain personal hygiene.</p>

100	Restoration of blood circulation	<p>A circulatory disorder is a condition in which blood does not flow in sufficient quantities to the brain, internal organs, and extremities. The problem can lead to severe complications, hypoxia, provoke deterioration of memory, hearing, etc. Such disorders develop amid some infectious or chronic diseases. The program helps to improve blood supply to organs and tissues especially after infectious diseases.</p> <p>Application: daily for a month, in conjunction with the programs on DeVita Ritm (see Heart and blood vessels complex).</p>
101	Healthy spleen	<p>The spleen is the organ responsible for hematopoiesis, immunity, and blood supply. Being in a constant state of stress, irritability leads to problems with the functioning of this organ, and usually leads to a decrease in immunity. The program supports spleen function by reducing the adverse effects of negative psychological factors .</p> <p>Application: daily for a month, in conjunction with DeVita Ritm programs (see "Stop stress" complex). In order to maintain immunity under constant stress, the scheme is used in a cycle of 2-3 times a year.</p>
102	Healing	<p>Healing of the body is impossible without the process of regeneration, which occurs at different levels of organization - system, organ, tissue, cellular, intracellular. The program activates the regeneration process in various inflammations (infectious, non-infectious) or damage to the integrity of organs and tissues (trauma), thereby enhancing the healing process.</p> <p>Application: daily for one month. If necessary, the program is used until full recovery.</p>
103	Stop Gardnerella	<p>Gardnerella is a bacterium that is part of the normal biocenosis of women, actively multiplying when local immunity is weakened, causing inflammation of the genitourinary organs. The program induces the death of bacteria, and reduces inflammation.</p> <p>Application: daily for 2-3 weeks, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see "Women's Health" complex). To prevent re-infection, sexual partners should use the program at the same time.</p>

104	Male power	<p>Men have faced the problem of erectile dysfunction at all times and regardless of age. One of the causes of this problem is chronic prostatitis. The program reduces the inflammatory process, contributes to the restoration of normal functioning of the reproductive system, reduces the chance of complications (impotence, infertility).</p> <p>Application: daily for a month, in conjunction with programs on DeVita Ritm (see "Men's Health" complex)"</p>
105	Without actinomycetes	<p>Actinomycetes are soil microorganisms. They form colonies on skin and mucous membranes without causing disease. They are capable of penetrating the human body only through damaged skin, mucous membranes. They become active under the influence of various factors, which include damage to the mucosa, the presence of infectious diseases, weakened immunity, diabetes mellitus, cancer, etc. The program prevents the spread of infection, and reduces the risks of complications.</p> <p>Application: daily, in conjunction with the "Without toxins", "Restoration 2- Automatic" programs on DeVita AP .</p>
106	Stop gonococcus	<p>A bacterium that causes inflammation of the genitourinary tract. It's transmitted sexually. The program inhibits microorganisms, and reduces inflammation.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see "Women's Health", "Men's Health" complexes). To prevent re-infection, sexual partners should use the above scheme at the same time.</p>
107	Stop parvovirus	<p>A virus that affects the joints and causes chronic anemia. The virus is dangerous for people with weakened immune systems and blood diseases. The program induces the death of the virus, and reduces the risk of complications.</p> <p>Application: daily month, in conjunction with the "Without toxins" program on DeVita AP.</p>
108	Without tuberculosis	<p>A chronic infectious disease caused by bacteria. It affects respiratory organs, bones and joints, genitourinary organs, eyes, peripheral lymph nodes. The following factors can provoke the disease: reduced immune defense, stress, chronic diseases of the lungs and GI tract, bad habits, living conditions. The program fights the root cause of the disease - the pathogenic flora.</p>

108	Without tuberculosis	Application: daily for a month, in conjunction with the "Without toxins" program on the DeVita AP and the "Healthy day" and "Health" programs on the DeVita Ritm.
109	Without brucella	Brucella is the bacteria that causes brucellosis. The disease is transmitted by ingestion of contaminated food (thermally unprocessed meat and dairy products), contact with an infected animal, or by aerosol transmission. The disease affects almost all organs, occurs in both acute and chronic forms. The program suppresses bacteria and reduces the toxic load on the body. Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm.
110	Skin without acne	Acne is a disease of the hair follicles and sebaceous glands that is more common in adolescents and young adults. There is a distinction between non-inflammatory and inflammatory acne. The program reduces inflammation, promotes the regulation of normal skin function, and improves the skin's condition. Application: daily for a month, in conjunction with the "Stop acne bacteria" program on DeVita AP and the corresponding programs on DeVita Ritm (see Clear Skin Complex). In case of skin predisposition to acne formation, in order to prevent inflammation, it is desirable to use the program every quarter according to the above scheme.
111	Fungus-free skin	Dermatomycosis is an infectious disease of the skin and its appendages, which is caused by various pathogenic, conditionally pathogenic fungi. The following factors contribute to the development of the disease: dermatitis, chronic diseases, poor personal hygiene, and professional factors. The program reduces toxic load and inflammation. Application: daily for a month, combined with anti-fungal programs.
112	Straight spine	Spinal curvature is a disruption in the proper alignment of the spinal column. Any curvature of the spine can lead to not only back pain, but also internal organ problems. Various factors influence the curvature of the spine: weakness of the muscular corset of the back; sedentary behavior; asymmetrical types of physical activity; constant tilting of the torso forward when sitting; carrying a bag or backpack on one side; diseases of the musculoskeletal system.

112	Straight spine	<p>The program helps reduce muscle spasm, reduces pain, promotes muscle recovery, and improves the health of all tissues of the spine.</p> <p>Application: daily, in conjunction with the "Healthy spine-automatic", "Healthy spine-automatic" programs on DeVita Ritm.</p>
113	Stop pain (low back)	<p>Discomfort in the lower back is common to almost every person in their 30s and older. There are many reasons for this: circulatory disorders; degenerative changes; trauma; infections, etc. The program reduces pain and improves blood flow in the lower back area.</p> <p>Application: daily for a month or as needed in conjunction with the "Low back without inflammation" program on the DeVita AP and the "Joints without pain-automatic" and "Healthy spine-automatic" programs on the DeVita Ritm.</p>
114	Without allergies new	<p>The occurrence of allergies in humans is influenced by a number of factors. Among them: household allergens, working conditions, genetic predisposition and other factors. These conditions can not only cause allergic reactions, but can contribute to the development of asthma. The program reduces allergies, promotes better breathing by relieving spasm, and it improves well-being.</p> <p>Application: daily for a month or as needed, in conjunction with the corresponding programs on program on DeVita AP and the DeVita Ritm (see Stop allergy complex).</p>
115	Without fluke	<p>The Siberian or cat liver fluke a flatworm which belongs to the class of trematodes. It causes opisthorchiasis, affecting organs such as the gallbladder, pancreatic ducts, and bile ducts of the liver. Infection occurs by eating low-salted, undercooked food, raw fish; by using household appliances that have come into contact with such fish; by failing to maintain personal hygiene after visiting a restroom or cleaning animal toilets and interacting with them. The program deals with the parasite at all stages of its development.</p> <p>Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the corresponding programs on DeVita Ritm (see Stop parasites complex). It is mandatory to maintain personal hygiene. To consolidate the result, it is desirable to repeat the cycle in a month.</p>

116	Stop dementia	<p>A number of reasons can lead to impaired thinking abilities in a person: vascular diseases, Alzheimer's disease, trauma, brain neoplasms, alcoholism, drug addiction, CNS infections and some other diseases. The program helps maintain brain function.</p> <p>Application: daily for a month, repeat the cycle at regular intervals every quarter. Use in conjunction with the "Blood circulation in the head" program and the corresponding programs on the DeVita Ritm.</p>
117	Stop Parkinson	<p>The disease is characterized by damage to nerve cells in the brain stem and is expressed by slowness of movement, tremors at rest and impaired reflexes. It is common in the elderly. The program helps to slow the progression of the disease.</p> <p>Application: daily for a month, repeat the cycle at regular intervals every quarter. Use in conjunction with the "Blood circulation in the head" program and the corresponding programs on the DeVita Ritm.</p>
118	Stop pain (injury)	<p>Nobody is safe from injury. Injuries are accompanied by pain, swelling, violation of tissue integrity and other signs. The program reduces pain, swelling, and speeds up the regeneration process.</p> <p>Application: daily for a month, in conjunction with the "Stop Injuries" and "Blood circulation regulation" program on DeVita Ritm.</p>
119	Strong nerves	<p>Amid a prolonged nervous overstress due to prolonged illness, life troubles, loss of loved ones and other reasons can arise nervous disorders characterized by irritability, general fatigue, apathy, mood swings. The program helps to normalize the psychological state, helps to improve mood and general condition.</p> <p>Application: daily for a month, in conjunction with DeVita Ritm programs (see "Stop Stress" complex).</p>
120	Good sleep	<p>Healthy sleep is the key to good mood and well-being. But stressful situations at home and work, lifestyle changes and other reasons can trigger insomnia. The program reduces the irritability of the nervous system and promotes restful sleep.</p> <p>Application: daily for a month, in conjunction with the "Strong nerves" program and the corresponding programs on DeVita Ritm (see Stop insomnia).</p>

121	Calm intestines	<p>Stomach heaviness, heartburn, bloating, rumbling, meteorism, stool problems can all be signs of irritable bowel. Such a condition is typical for people in developed countries. The program reduces the severity of pain, relieves mucosal inflammation, normalizes GI motility.</p> <p>Application: daily for a month, in conjunction with the "Intestines without disorders" program on DeVita AP and the corresponding programs on DeVita Ritm (see Digestion without problems complex).</p>
122	Healthy Stomach	<p>Stomach pain is familiar to almost every other person. The causes of such pains are: functional disorders, intestinal infections, gastritis, etc. Stress, improper nutrition, taking some medications can provoke pain. The program relieves pain, reduces inflammation, helps normalize the stomach.</p> <p>Application: daily for a month, in conjunction with the "Healthy stomach" program on DeVita Ritm and the "Stop stress" or "Digestion without problems" complex programs. In order to prevent exacerbations, the cycle is carried out 2-3 times a year.</p>
123	Duodenitis, stop!	<p>Eating disorders, consumption of foods that irritate the mucosa of the GI tract, diseases of the stomach, liver and pancreas can cause inflammation of the mucous membrane of the duodenum and pancreas, which is often accompanied by severe pain. The program relieves pain and helps mucous membrane recovery.</p> <p>Application: daily for a month, in conjunction with the corresponding programs on the DeVita AP and the "Stomach+ Automatic" program on the DeVita Ritm.</p>
124	Healthy pancreas	<p>The pancreas plays an important role in the digestive and endocrine systems, it produces enzymes that aid the digestive process and the hormone insulin that regulates blood sugar levels. Therefore, it is very important to keep it functioning properly. The program reduces discomfort and maintains the function of the pancreas.</p> <p>Application: daily for a month, in conjunction with the corresponding programs on DeVita AP and DeVita Ritm.</p>
125	Blood sugar under control	<p>Stress, high-carbohydrate diets, medications, and endocrine and internal organ diseases lead to high blood sugar. It is possible to balance sugar levels by combining diet, exercise and the "Blood sugar under control" program.</p> <p>Application: daily for one month.</p>

126	Stop diabetes	<p>Diabetes is a disease characterized by chronically elevated blood glucose levels. There is a risk of cardiovascular disease, loss of vision, decreased immunity, and impaired blood supply to the lower extremities. The program helps to minimize possible risks.</p> <p>Application: daily for a month, in conjunction with the program "Blood sugar under control" on DeVita AP and basic programs on DeVita Ritm complex "Digestion without problems" and special programs "Normal blood sugar", "Healthy Pancreas". The full cycle is performed 2 times a year, the rest of the time the DeVita Ritm programs are used.</p>
127	Small intestine	<p>The small intestine performs several important functions: digestion of food; absorption of nutrients; movement of remaining substances to the GI tract; hormone production; and immunological defense. Most often people are faced with functional disorders and inflammatory processes in this part of the intestine. The program helps to reduce pain and inflammatory manifestations, helps to normalize bowel function.</p> <p>Application: daily for a month, in conjunction with the corresponding programs on program on DeVita AP and the DeVita Ritm (see Digestion without problems complex).</p>
128	Intestines without disorders	<p>Disorders of the intestines can appear with functional impairments in its work, infectious diseases or stress. They are accompanied by abdominal pain, bloating, and constipation and/or diarrhea. The program reduces pain, helps normalize the work of the large intestine.</p> <p>Application: daily for a month, in conjunction with the "Calm intestines" program on DeVita AP and the corresponding programs on DeVita Ritm (see Digestion without problems complex).</p>
129	Without hemorrhoids	<p>Hemorrhoids are the result of impaired blood circulation and inflammation in the blood vessels of the rectum. The disease is equally common in men and women, and brings great discomfort. The program relieves pain, reduces inflammation, and improves vascular health.</p> <p>Application: daily for one month. To minimize the occurrence of the disease, it is necessary to control the influence of risk factors; take a cycle using the schemes of complexes "Digestion without problems", "Heart and blood vessels" or others, depending on the cause of the problem.</p>

130	Stop pelvic inflammation	<p>Inflammatory processes in the pelvic area can be both independent diseases and complications of previously experienced diseases. One of the most common causes can be considered bacterial and viral infections. The program reduces inflammation and relieves pain.</p> <p>Application: daily for a month, in conjunction with the DeVita AP and DeVita Ritm programs for relevant infectious agents (see "Women's Health" and "Men's Health").</p>
131	Stop chickenpox	<p>Chickenpox is an infectious disease caused by the varicella-zoster virus that affects the skin and nervous system. It is transmitted by airborne and contact. It is more severe in adults. The program fights the pathogen, reducing the manifestation of intoxication and the occurrence of complications.</p> <p>Application: daily for 2 weeks, in conjunction with the "Without toxins" program on DeVita AP.</p>
132	Stop rubella	<p>Rubella is an acute viral infection that is very "contagious". It's transmitted by airborne droplets. The program inhibits the development of the virus and its spread in the body.</p> <p>Application: daily month, in conjunction with the "Without toxins" program on DeVita AP.</p>
133	Without mumps	<p>Mumps is an acute infectious disease caused by a virus of the paramyxovirus genus. The main signs of the disease are intoxication and inflammation of the parotid salivary glands. The virus is transmitted by airborne and household contact. The program fights the virus, reducing the manifestation of intoxication and the occurrence of complications.</p> <p>Application: daily for 2 weeks, in conjunction with the "Without toxins" program on DeVita AP.</p>
134	Stop scarlet fever	<p>Scarlatina is an acute bacterial infection that is transmitted by airborne and household contact. Characterized by disease of the upper respiratory tract, tonsils. The program destroys pathogenic bacterial flora, reducing intoxication, reducing the risks of complications.</p> <p>Application: daily month, in conjunction with the "Without toxins" program on DeVita AP.</p>
135	Stop measles	<p>Rubella is an acute viral infection that is very "contagious". It's transmitted by airborne droplets. The program suppresses the development of the virus, relieves signs of intoxication, reduces the chances of complications.</p>

135	Stop measles	Application: daily month, in conjunction with the "Without toxins" program on DeVita AP.
136	No diphtheria	Diphtheria is an acute bacterial infection that is transmitted by airborne droplets and airborne dust. It is possible to be infected through the contact and domestic means (utensils, dirty hands, etc.), as well as alimentary (food). The disease is characterized by lesions of the upper respiratory tract, and oropharyngeal mucosa. The program fights the pathogen, reduces signs of intoxication and prevents the spread of inflammation to other organs. Application: daily month, in conjunction with the "Without toxins" program on DeVita AP. Use the appropriate programs on the DeVita Ritm and DeVita AP if necessary.
137	Stop cholera	Cholera is an acute bacterial infection, transmitted by fecal-oral route through water or contaminated food. The main symptoms are vomiting, diarrhea, and abdominal pain. The program inhibits the development of the bacteria, which helps to reduce intoxication and improve the person's condition. Application: daily month, in conjunction with the "Without toxins" program on DeVita AP. It is mandatory to maintain personal hygiene. In order to prevent the disease, for persons visiting countries in Asia (India), Africa, South America, each time during and/or after travel - daily for 2 weeks.
138	Stop pertussis	Pertussis is an acute bacterial infection with damage to the respiratory, cardiovascular and nervous systems. The path of transmission of the pathogen is airborne. The program fights the pathogen by reducing the effects of toxins on cells and preventing inflammation from spreading to other organs. Application: daily for a month, in conjunction with the "Without toxins" program on DeVita AP and the "Cellular respiration" and "Free Breath" programs on DeVita Ritm.
139	Without constipation	Improper nutrition (lack of dietary fiber, liquids), sedentary work, stress, change of climatic area during relocation, impaired motility of the large intestine and other factors can cause difficult or systematically obstructed defecation. Prolonged chronic constipation can contribute to the development of complications in the digestive tract. The program helps to stimulate bowel function and strengthen the intestinal muscle tone.

139	Without constipation	<p>Application: daily for a month, in conjunction with the corresponding programs on DeVita AP and DeVita Ritm (see "Digestion without problems", "Stop stress" complexes. To prevent further constipation, introduce fiber-rich foods into your diet and eliminate foods that aggravate constipation. Increase fluid intake. With a predisposition to this problem, frequent relocations, stress, the cycle is carried out 2-3 times a year.</p>
140	Blood circulation in the head	<p>Insufficient blood supply to the organs of the head can lead to undesirable consequences: headache, dizziness, tinnitus, decreased visual and hearing acuity, impaired perception, attention, memory. The program helps to improve blood flow to the tissues of the head organs, thereby increase the flow of oxygen and nutrition to them.</p> <p>Application: daily for a month, in conjunction with the corresponding programs on program on DeVita AP and the programs DeVita Ritm (see Heart and Blood Vessels complex).</p>
141	For women	<p>Abundant and painful menstruation, irregularity of the cycle, may be the cause of endometriosis. One of the complications of this disease is infertility. There are many causes of the disease: hormonal disorders, decreased immunity, heredity, etc. The program helps to restore hormonal balance, normal functioning of the reproductive system.</p> <p>Application: daily for a month, in conjunction with the "Stop pelvic inflammation" program on DeVita AP and the programs DeVita Ritm (see Women's health complex).</p>
142	Stop spasm (muscles)	<p>The superficial and deep back muscles are like a corset, supporting the spine in the correct position. Lifting heavy weights, spending long periods of time in an uncomfortable position, poor posture, problems with the musculoskeletal system cause spasms of the back muscles, accompanied by pain. The program helps to relieve spasm, pain, reduces inflammation, and improves a person's condition.</p> <p>Application: daily for a month, in conjunction with the DeVita AP "Muscles without tension" program and the DeVita Ritm "Healthy spine-automatic" and "Joints without pain-automatic" programs. To strengthen the muscular corset of the back, regular training of all back and abdominal muscles using special DeVita Ritm programs is essential (see "For Athletes" complex).</p>

143	Without motion sickness	<p>Motion sickness is caused by the vestibular apparatus – the organ of balance and orientation in space, which is located in the inner ear. When motion sickness occurs, feelings of nausea and dizziness are experienced. The program relieves cramps, nausea, and improves health.</p> <p>Application: every time you travel in transport, in conjunction with the “Blood circulation regulation” program on the DeVita Ritm.</p>
144	Stop spasm (Intestines)	<p>The pain associated with intestinal spasm can occur in connection with violations of various mechanisms of intestinal motility, with intoxication, against the background of digestive disorders, stressful situations, helminthiasis, etc. The program helps to relax the intestinal wall, reduce pain.</p> <p>Application: daily for a for one month, in conjunction with the corresponding programs on DeVita AP and DeVita Ritm (see “Digestion without problems”, “Stop stress” complexes).</p>
145	Stop conjunctivitis	<p>Conjunctivitis is an inflammation of the mucous membrane, occurs most often with viral or bacterial infection, allergies. It is characterized by redness, swelling and itching. The program helps to reduce inflammation, swelling and itching.</p> <p>Application: for infections daily for a month, in conjunction with the appropriate programs on the DeVita AP . For allergies, daily for the entire period of exposure to allergens, in conjunction with the corresponding programs on DeVita AP and DeVita Ritm (see “Stop Allergy” complex).</p>
146	Muscles without tension	<p>Muscle tension can be caused by injuries, uncomfortable posture, poor posture, osteochondrosis, stress, or disease. The program reduces pain, increases muscle elasticity, and reduces muscle tone.</p> <p>Application: daily for a month, in conjunction with the “Stop spasm (muscles) program on DeVita AP and “Muscles without pain” on DeVita Ritm.</p>
147	Stop fatigue	<p>Fatigue is characterized as a state of exhaustion, manifested by a feeling of weakness, lethargy, feelings of physiological discomfort, disturbances in the flow of mental processes. Fatigue can be the result of overwork, poor sleep, worries, lack of exercise, anxiety, and illness. It is most frequent in residents of large cities. The program reduces the irritability of the nervous system and helps to increase performance.</p>

147	Stop fatigue	Application: daily for a month, in conjunction with the "Strong nerves", "Blood circulation in the head" programs on DeVita AP and the programs DeVita Ritm (see Stop stress complex).
-----	--------------	---

www.deholding.info

Deta-Elis Europa GmbH
Kaiserstrasse 74, 63065 Offenbach, Germany